




petrotest[®]

PETROLEUM TEST EQUIPMENT PRODUCT RANGE 2010 / 2011



MANUFACTURER SINCE 1873

WWW.PETROTEST.COM



Contact

Mail:

Petrotest@ Instruments
GmbH & Co. KG
Ludwig-Erhard-Ring 13
15827 Dahlewitz
Germany

e-Mail:

info@petrotest.com

Website:

www.petrotest.com

Fax:

+49 (0)33 708 / 56 -556

Phone:

+49 (0)33 708 / 56 -300

Direct line to:

- Experts: -350
- Service: -442
- Sales: -521
- Shipping: -521
- Marketing: -348
- Purchase: -336
- Accounting: -344
- Management: -314

ISO 9001 certified

Company Development

From Berthold Pensky

In 1873, Mr. Berthold Pensky, an engineer at the "Kaiser-Wilhelm-Institute" in Berlin-Dahlem, today known as the "Bundesanstalt für Materialprüfung" (German Federal Institution for Materials Testing), founded his own company and produced Pensky-Martens flash-point testers, Abel-Pensky flash-point testers, as well as viscometers, grain testers and other machinery.



1909-1989

..... through Sommer & Runge

1886, upon his appointment to the official position of a "Royal Prussian Planning Officer" Mr. Pensky sold his company to two of his employees, Mr. Wilhelm Sommer and Mr. Max Runge. By then, the above products were marketed in Germany, the Netherlands, Russia and Romania already. During the bombing of Berlin the factory was partially destroyed and had to move to other quarters. In 1948 a manufacturing line of mineral oil test instruments for the German speaking market was started.

..... to Petrotest® Instruments

Following the reunification of Germany, Sommer & Runge (SuR) reformed and renamed under Petrotest® Instruments GmbH & Co KG. Due to the need for expansion Petrotest® relocated in a brand new building on the outskirts of Berlin in the Dahlewitz Industrial Park.

Approximately 60 employees are engaged in the production of mechanical and electronic laboratory instruments such as flash-point testers, Fraass breaking point testers, penetrometers etc. Further improvement and development of more automatic petroleum testing equipment can be expected.



since 1996

Petrotest® Philosophy



In order to be a highly flexible and reliable partner to our customers, we decided to locate the entire range of branches needed in our company.

- Application
- Customer Care Centre
- Electrical & Mechanical Design
- Software Development
- Production incl. CNC-machinery
- Assembling
- Quality Control according ISO 9001
- Marketing
- Order & Shipping
- Stockage

This gives us the potentiality to provide even customized solutions. Therefore do not hesitate to ask.



Chapters

	Page
(A)  Sampling Devices	4
(B)  Density & Gravity Testing	5 - 6
(C)  Viscosity Testing	7 - 9
(D)  Distillation Units	10
(E)  Flash-Point & Auto-Ignition Testing	11 - 13
(F)  LPG Testing	14
(G)  Fuel Testing	15 - 19
(H)  Constituent Analyzer (S / Cl / Salt etc.)	20 - 22
(K)  Oil Testing	23 - 30
(K)  Field Tests - Quick Checks	30
(L)  Friction & Wear Testing	31
(M)  Grease & Wax Testing	32 - 35
(N)  Bitumen Testing	36 - 38
(P)  Thermometers	39
(Q)  Calibration Fluids	39
(R)  Metal Specimen	39
(T)  Penetration & Texture Testing	40
(V)  Various Equipment	41 - 42



ADU



PMA



BPA



PNR

Tables of Contents

by Title	44 - 47
by Order Numbers	48 - 49
by Standards	50 - 53



Sampling Devices



19-0051

ASTM D 923 - ASTM D 4057 - ASTM D 4177 - DIN 51 750-1-3 - EN 58 - FTM 791-8001 - ISO 3170 - ISO 15 528 - JIS C2101 - JIS K2254 (ASTM D 270, DIN 53 242-1-4 obs.)

Datasheet:
19-.pdf



19-0056

Sampling Tubes - open (Grease)

ASTM D 4057 (ASTM D 270 obs.) - DIN 51 750-3 - EN 58 (DIN 53 242-3, IP 51 obs.)

Ref.-No:
19-0051 (Auger)
19-0056 (ASTM)
19-0058 (DIN)



19-0058



19-0005

Sampling Tubes - open

DIN 51 750-2 - EN 58 - ISO 3170

Ref.-No:
19-0005



19-0002

Sampling Tubes - closed

DIN 51 750-1+2 - EN 58 - ISO 3170

Ref.-No:
19-0001
19-0002
19-0003
19-0011
19-0012



19-0032 19-0040 19-0041

Tank Sampling Bottles

ASTM D 923 - ASTM D 4057 (ASTM D 270 obs.) - DIN 51 750-2 - EN 58 - ISO 3170 - (DIN 53 242-1 obs.) - JIS C2101 - JIS K2254

Ref.-No:
19-0032 (BACON)
19-0040 (Can)
19-0041 (Glas)



Tank Sampling by Scoop-Thermometer

scoop Ø 35 mm, length 100 mm
19-0070 with thermometer range 0...+60°C
19-0071 with thermometer range 0...+100°C

Datasheet:
19-.pdf

Ref.-No:
19-0070
19-0071



Lowering Equipment

Wheel with Wire Cord (19-0080)
Brass Chain (19-0035 to 19-0037)
Hemp Rope (19-0034)
Perlon Cord (19-0039)



19-0153 & 19-0154

Sampling Cylinders - double valve

ASTM D 1265
Volume: 75 ml to 3785 ml
acc. DOT 3E or 3A requirements

Datasheet:
19-0150.pdf

Ref.-No:
19-0150
19-0151
19-0152
19-0153
19-0154
19-0155
19-0156



Tempering Bath for Density, API & Specific Gravity

ASTM D 70 - ASTM D 71 - ASTM D 287 - ASTM D 1298 -
 ASTM D 1481 - ASTM D 2111 - ASTM E 100 - DIN 51 757 - DIN 52 004
 (DIN 1995 obs.) - IP 160 - IP 189 - IP 190 - ISO 3675 - ISO 3838 -
 JIS K2207 - JIS K2249 - JIS K2265 (BS 4714 - BS 4699 obs.)

- 4-place full-view acryl bath
- for determination of Density, Specific & API Gravity
- up to +60 °C max.

Datasheet:
 26-0010.pdf

Ref.-No:
26-0010 (230V)
26-0011 (115V)



Tempering Bath for Density, API & Specific Gravity

ASTM D 70 - ASTM D 71 - ASTM D 287 - ASTM D 1298 -
 ASTM D 1481 - ASTM D 2111 - ASTM E 100 - DIN 51 757 - DIN 52 004
 (DIN 1995 obs.) - IP 160 - IP 189 - IP 190 - ISO 3675 - ISO 3838 -
 JIS K2207 - JIS K2249 - JIS K2265 (BS 4714 - BS 4699 obs.)

- 4-place full-view borosilicate bath
- for determination of Density, Specific & API Gravity
- up to +150 °C max.

Datasheet:
 26-0017.pdf

Ref.-No:
26-0017 (230V)
26-0018 (115V)



Hydrometer, DIN 12 791 - Series M 100

DIN 51 757

Ref. temp.: +20 °C, max. length: 250 mm, scale division: 0.0020 g/cm³

Datasheet:
 26-0400 etc.pdf

Ref.-No:
26-0400
 to
26-0405

Hydrometer, API Specific Gravity - ASTM E 100

ASTM D 287 - ASTM D 1298 - ASTM E 100 - IP 160 - ISO 3675 -
 JIS K2249

Ref. temp.: +60 °F, max. length: 330 mm, scale division: 0.1 API

Datasheet:
 26-0400 etc.pdf

Ref.-No:
26-0500
 to
26-0509



Hydrometer, ASTM Specific Gravity - ASTM E 100

ASTM D 1298 - ASTM E 100 - IP 160 - JIS K2249

Ref. temp.: +60 °F, max. length: 330 mm, scale division: 0.0005 g/cm³

Datasheet:
 26-0400 etc.pdf

Ref.-No:
26-0600
 to
26-0608

Hydrometer, Density - Series M 50 SP

ASTM D 1298 - IP 160 - JIS K2249

Ref. temp.: +15 °C, max. length: 270 mm, scale division: 0.001 g/cm³

Datasheet:
 26-0400 etc.pdf

Ref.-No:
99-1141
 to
99-1150



Hydrometer, Density - Series L 50 SP

ASTM D 1298 - IP 160 - JIS K2249

Ref. temp.: +15 °C, max. length: 335 mm, scale division: 0.0005 g/cm³

Datasheet:
 26-0400 etc.pdf

Ref.-No:
99-1131
 to
99-1140

Hydrometer, Density - Series L 50 SP

ASTM D 1298 - IP 160 - JIS K2249

Ref. temp.: +20 °C, max. length: 335 mm, scale division: 0.0005 g/cm³

Datasheet:
 26-0400 etc.pdf

Ref.-No:
99-1031
 to
99-1040



Hydrometer, Specific Gravity - Series M 50 SP

ASTM D 1298 - IP 160 - JIS K2249

Ref. temp.: +60 °F, max. length: 270 mm, scale division: 0.001 g/cm³

Datasheet:
26-0400 etc.pdf

Ref.-No:
99-1151
to
99-1160



Hydrometer, Specific Gravity - Series L 50 SP

ASTM D 1298 - IP 160 - JIS K2249

Ref. temp.: +60 °F, max. length: 335 mm, scale division: 0.0005 g/cm³

Datasheet:
26-0400 etc.pdf

Ref.-No:
99-1161
to
99-1170



Hydrometer (Thermo-), Specific Gravity - ASTM 101H

ASTM D 1657

Datasheet:
26-0400 etc.pdf

Ref.-No:
14-0116 (°F)
14-0121 (°C)



Pycnometer - LIPKIN

ASTM D 1481 - IP 189 - ISO 3838 (DIS 3658 obs.) -
JIS K2249 (BS 4699 obs.)

Datasheet:
26-0712 etc.pdf

Ref.-No:
26-0712



Pycnometer - GAY-LUSSAC

IP 190 (BS 4699 obs.) - ISO 3838 (DIS 3658 obs.) - JIS K2249

Datasheet:
26-0712 etc.pdf

Ref.-No:
26-0721



Pycnometer - HUBBARD

ASTM D 70 - ISO 3838 (DIS 3658 obs.) - JIS K2249

ASTM D 70 - DIN 52 004 (DIN 1995 obs.) - ISO 3838 (DIS 3658 obs.) -
IP 190 (BS 4699 obs.) - JIS K2249

Datasheet:
26-0712 etc.pdf

Ref.-No:
26-0740 (IP)
26-0741 (ISO)

page 14

Density or Relative Density of Light Hydrocarbons by Pressure Hydrometer

see: LPG

page 14

Density or Relative Density of Gases by SCHILLING Effusimeter

see: LPG



Density Balance (digital)

ASTM D 71 - DIN 51 757-2 (IP 59-A+B obs.)

- conversion of measured balance to density
- build-in RS 232 C interface
- platform with solid stand

Datasheet:
26-0019.pdf

Ref.-No:
26-0019 (115/230V)



**Tempering Baths for Viscosity
- Kinematic Capillary Methods**

ASTM D 445 - ASTM D 446 - ASTM D 2170 - (DIN 51 561, DIN 51 569 obs) - EN 12 595 - FTM 791-305 - IP 71-1+2 - IP 319 - ISO 3104 (NF T60-100 obs.) - ISO 3105 - JIS K2207 - JIS K2283

**Tempering Bath for Viscositytubes
- Ambient Temperature**

- 4-place full-view acryl bath
- max. temperature up to +60 °C
- for Ubbelohde-, Cannon-Fenske- or BS/IP-U-Tube-viscometers

Datasheet:
20-0000.pdf

Ref.-No:
20-0000 (230V)
20-0001 (115V)



**Tempering Bath for Viscositytubes
- High Temperature**

- 7-place bath
- temperature up to +150 °C max.
- for Ubbelohde-, Cannon-Fenske- or BS/IP-U-Tube-viscometers

Datasheet:
20-0150.pdf

Ref.-No:
20-0150 (230V)
20-0151 (115V)



Viscosity by UBBELOHDE Capillary - AKV (automatic)

ASTM D 445 - ASTM D 446 - IP 71-1+2 - JIS K2283

Datasheet:
22-1000.pdf

Ref.-No:
22-1000 (230V)
22-1001 (115V)



**Viscosity by UBBELOHDE ASTM/DIN Capillary
- ViscoClock (semi-automatic)**

ASTM D 445 - ASTM D 446 - DIN 51 562 - IP 71-1+2 - IP 319 - ISO 3104 (NF T60-100 obs.) - ISO 3105 - JIS K2283

- fits almost every customer's tempering bath
- for bath temperature ranges of -40 to +150 °C
- interface for PC or printer

Datasheet:
22-0175.pdf

Ref.-No:
22-0175 (230V)
22-0176 (115V)



**Viscometer, UBBELOHDE ASTM-Type
- UAV (Suspended Level)**

ASTM D 445 - ASTM D 446 - FTM 791-305 - IP 71-1+2 - ISO 3104 (NF T60-100 obs.) - ISO 3105 - JIS K2283

Datasheet:
22-0120.pdf

Ref.-No:
22-0120 (manual)
22-0320 (automatic)



**Viscometer, UBBELOHDE DIN-Type
- UDV (Suspended Level)**

ASTM D 445 - ASTM D 446 - DIN 51 562-1 - DIN 52 007-1 - ISO 3104 (NF T60-100 obs.) - ISO 3105 - JIS K2283

Datasheet:
22-0100.pdf

Ref.-No:
22-0100 (manual)
22-0300 (automatic)



**Viscometer, UBBELOHDE Micro Viscometers
- UMV (Suspended Level)**

ASTM D 445 - ASTM D 446 - DIN 51 562-2 - JIS K2283

Datasheet:
22-0390.pdf

Ref.-No:
22-0390 (manual)
22-0395 (automatic)



**Viscometer, CANNON-UBBELOHDE Semi-Micro
- CUSMV (Suspended Level)**

ASTM D 445 - ASTM D 446 - IP 71-1+2 - ISO 3104 (NF T60-100 obs.) - ISO 3105 - JIS K2283

Datasheet:
22-0360.pdf

Ref.-No:
22-0360 (manual)





Viscometer, CANNON-FENSKE Routine - CFRV (Ostwald)

ASTM D 445 - ASTM D 446 - ASTM D 2170 - FTM 791-305 -
IP 71-1+2 - IP 319 - ISO 3104 (NF T60-100 obs.) - ISO 3105 -
JIS K2283

Datasheet:
22-0000.pdf

Ref.-No:
22-0000 (manual)



Viscometer, CANNON-FENSKE Opaque - CFOV (Reverse Flow)

ASTM D 445 - ASTM D 446 - ASTM D 2170 - DIN 51 366 -
(DIN 51 372 obs., DIN 51 550 obs.) - DIN 52 007-1 - EN 12 595 -
FTM 791-305 - IP 71-1+2 - IP 319 - ISO 3104 (NF T60-100 obs.) -
ISO 3105 - JIS K2207 - JIS K2283

Datasheet:
22-0050.pdf

Ref.-No:
22-0050 (manual)



Viscometer, BS/IP U-Tube - BS/IP/RF (Reverse Flow)

ASTM D 445 - ASTM D 446 - ASTM D 2170 - DIN 51 366 -
(DIN 51 372 obs., DIN 51 550 obs.) - DIN 52 007-1 - EN 12 595 -
FTM 791-305 - IP 71-1+2 - IP 319 - ISO 3104 (NF T60-100 obs.) -
ISO 3105 - JIS K2207 - JIS K2283

Datasheet:
22-0200.pdf

Ref.-No:
22-0200 (manual)



Viscometer, BS/IP - BS/IP/SL(S)(Suspended Level)

ASTM D 445 - ASTM D 446 - FTM 791-305 - IP 71-1+2 -
ISO 3104 (NF T60-100 obs.) - ISO 3105 - JIS K2283

Datasheet:
22-0230.pdf

Ref.-No:
22-0230 (manual)



Viscometer, Micro-OSTWALD - MOV (Ostwald)

ASTM D 445 - ASTM D 446 - FTM 791-305 - IP 71-1+2 -
ISO 3104 (NF T60-100 obs.) - ISO 3105 - JIS K2283

Datasheet:
22-0400.pdf

Ref.-No:
22-0400 (manual)
22-0405 (automatic)



Viscometer, ZEITFUCHS Cross-Arm - ZCAV (Reverse Flow)

ASTM D 445 - ASTM D 446 - ASTM D 2170 - EN 12 595 -
FTM 791-305 - IP 71-1+2 - IP 319 - ISO 3104 (NF T60-100 obs.) -
ISO 3105 - JIS K2207 - JIS K2283

Datasheet:
22-0530.pdf

Ref.-No:
22-0530 (manual)

Viscosity, Oil Standards

ASTM D 2162 - ASTM D 4052 - ISO 9001 (DIN 51 550 obs.)

Datasheet:
24-0551.pdf

Ref.-No:
24-0551 (Primary)
24-0501 (Secondary)
24-0601



Viscometer, Vacuum - CMVV, AIVV, MKVV

ASTM D 2171 - EN 12 596 - IP 222

22-0500 Cannon-Manning - CMVV
22-0520 Viscometer, Asphalt Institute - AIVV
Viscometer, Modified Koppers - MKVV

Datasheet:
22-0500.pdf

Ref.-No:
22-0500
22-0520



Viscosity by REDWOOD (manual)

IP 70 (obs.)

21-0300 with oil cup n° 1
21-0335 with oil cup n° 2

Datasheet:
21-0300.pdf

Ref.-No:
21-0300 (230V)
21-0335 (230V)



Viscosity by SAYBOLT (manual)

ASTM D 88 - ASTM E 102 - AASHTO T72 - FTM 791-304 - JIS K2207

- 2 viscosity tubes, gold plated
- stainless steel housing
- thermostatically controlled test bath with electric stirrer

Datasheet:
21-0501.pdf

Ref.-No:
21-0501 (230V)
21-0502 (115V)



Viscosity by Falling Ball - HÖPPLER (manual)

DIN 53 015

Datasheet:
24-0025.pdf

Ref.-No:
24-0025



Viscosity by Falling Ball - JUNG Check

see: Oil

page 30

Viscosity by Comparator - Viscosity Check

see: Oil

page 30

Viscosity by Flow Cup - FORD (manual)

DIN 53 211

ASTM D 1200 (NF T30-070 obs.)

Datasheet:
21-0761.pdf

Ref.-No:
21-0750 (DIN)
21-0761 (ASTM)



Viscosity by Flow Cup - ISO (manual)

ASTM D 5125 - DIN 53 224 - ISO 2431 (EN 535 obs.)

Datasheet:
21-0761.pdf

Ref.-No:
21-0766 (ISO)





Distillation at Atmospheric Pressure - ADU 4+ (automatic)

ASTM D 86 - IP 123 - ISO 3405 - JIS K2254
ASTM D 850 - ASTM D 1078 - IP 195

ADU 4+ : basic unit
ADU 4+ DryPoint : basic unit with dry-point detection

- built-in cooling system
- built-in fire-extinguisher
- integrated computer & external TFT monitor, keypad and mouse
- software included

Datasheet:
11-5580.pdf

Ref.-No:
11-5580 (230V)
11-5581 (115V)

11-5582 (230V)
11-5583 (115V)
- Dry-Point



Distillation at Atmospheric Pressure - ADU 4 (automatic)

ASTM D 86 - IP 123 - ISO 3405 - JIS K2254
ASTM D 850 - ASTM D 1078 - IP 195

ADU 4 : basic unit
ADU 4-DryPoint : basic unit with dry-point detection

- built-in fire-extinguisher
- software included
- ultra-high heating power
- safety-heater-rods reduce burn-risk

Datasheet:
11-0542.pdf

Ref.-No:
11-0542 (230V)
11-0543 (115V)

11-1542 (230V)
11-1543 (115V)
- Dry-Point



Distillation at Atmospheric Pressure - DU 4-Economy (manual)

ASTM D 86 - IP 123 - ISO 3405 - JIS K2254
ASTM D 850 - ASTM D 1078 - IP 195

DU 4-Eco
basic unit with: - 800W heater
- external bath-tempering
- for distillation groups 0 - 3 & 4 (up to +300°C)

DU 4-EcoPlus
basic unit with: - 1200W heater (safety-heater-rods reduce burn-risk)
- external bath-tempering
- for distillation groups 0 - 4

DU 4-EcoThermo
basic unit with: - 800W heater
- with built-in bath-heater but external cooling
- for distillation groups 0 - 3 & 4 (up to +300°C)

DU 4-EcoThermoPlus
basic unit with: - 1200W heater (safety-heater-rods reduce burn-risk)
- with built-in bath-heater but external cooling
- for distillation groups 0 - 4

Datasheet:
11-0900.pdf

Ref.-No:
11-0900 (230V)
11-0901 (115V)
- Eco

11-0904 (230V)
11-0905 (115V)
- EcoPlus

11-0902 (230V)
11-0903 (115V)
- EcoThermo

11-0906 (230V)
11-0907 (115V)
- EcoThermoPlus



Distillation at Atmospheric Pressure - DU 4-Professional (manual)

ASTM D 86 - IP 123 - ISO 3405 - JIS K2254
ASTM D 850 - ASTM D 1078 - IP 195

DU 4-Pro : basic unit with built-in fire-extinguisher
DU 4-ProThermo : .. with built-in fire-extinguisher & bath heating

- low price - high value
- available for distillation groups 0 - 3 or 0 - 4
- ultra-high 1200 W heater (safety-heater-rods reduce burn-risk)

Datasheet:
11-0900.pdf

Ref.-No:
11-0910 (230V)
11-0911 (115V)
- Pro

11-0912 (230V)
11-0913 (115V)
- ProThermo



Distillation of Pitch

ASTM D 2569

Datasheet:
10-0179.pdf

Ref.-No:
10-0179 (230V)



Flash-Point Testing with Closed Cup acc. ABEL

- ABA 4 (automatic)

- IP 170 - ISO 13 736
 NF T66-009
 ISO 3679 - ISO 3680
 DIN 51 755-1-2 - NF M07-036
 ISO 1516 - ISO 1523 - EN 924 - IP 491 - IP 492
 ABA 4 - Air : air cooled version (12-0501)
 ABA 4 - Water : water cooled, low temperature version (12-0503)
- 1-twist handling
 - ABAcon PC-software included
 - RAPID-flash set and milli-cup available

- Datasheet:
 12-0501.pdf
- Ref.-No:
12-0501 (115/230V)
 - air cooled
12-0503 (115/230V)
 - water cooled



Flash-Point Testing with Closed Cup acc. TAG

- TAG 4 (automatic)

- ASTM D 56 - FTM 791-1101 - JIS K2265
 ASTM D 3934 - ASTM D 3941 - EN 924
 ISO 3679 - ISO 3680
 ASTM D 3934 - ASTM D 3941 - EN 924 - IP 491 - IP 492 -
 ISO 1516 - ISO 1523
 TAG 4 - Air : air cooled version (12-0540)
 TAG 4 - Water : water cooled, low temperature version (12-0542)
- 1-twist handling
 - TAGcon PC-software included
 - milli-cup & RAPID accessories available

- Datasheet:
 12-0540.pdf
- Ref.-No:
12-0540 (115/230V)
 - air cooled
12-0542 (115/230V)
 - water cooled



Flash-Point Testing with Closed Cup acc. PENSKY-MARTENS - PMA 4 (automatic)

- ASTM D 93-A+B - JIS K2265 - AASHTO T73 - AASHTO T172 -
 FTM 791-1102 - FTM 141-4293 - IP 34-A+B - ISO 2719-A+B
 (EN 22 719 obs.)
 ASTM D 6751
 ISO 15 267 - BS 684-1.17
- 1-twist handling
 - PMAcon PC-software included
 - fire-extinguisher option

- Datasheet:
 12-1770.pdf
- Ref.-No:
12-1770 (115/230V)
12-1767 (115/230V)
 - with fire-extinguisher



Flash-Point Testing - Sample Changer with PMA 4 (automatic)

- PMA 4 S: standard version (12-2770 or 12-2771)
 PMA 4 FS: fire-extinguishing version (12-2772 or 12-2773)
 PMA 4 SC: chilled sampler set version (12-2776 or 12-2777)
 PMA 4 FSC: fire-extinguishing & chilled sampler set version
 (12-2778 or 12-2779)
- up to 12 test-places
 - software included
 - milli-cup option

- Datasheet:
 12-2770.pdf
- Ref.-No:
12-2770 (with EU-plug)
12-2771 (with US-plug)
12-2772 (with EU-plug)
12-2773 (with US-plug)
12-2776 (with EU-plug)
12-2777 (with US-plug)
12-2777 (with EU-plug)
12-2778 (with US-plug)



Flash- & Fire-Point Testing with Open Cup acc. CLEVELAND Open Cup - CLA 5 (automatic)

- ASTM D 92 - FTM 791-1103 - IP 36 - ISO 2592 - JIS K2265 -
 AASHTO T48
 Superpave®
- barometric pressure correction
 - skimmer & draft deflector available
 - bitumen test option
 - PC-software included

- Datasheet:
 12-5700.pdf
- Ref.-No:
12-5700 (110/230V)





Auto-Ignition Temperature - ZPA 3 (automatic)

EN 14 522 - DIN 51 794 - CEI 79-4 - NF T20-037
NF T20-036 - 79/831/EWG

- autosampler available with or without heated sample-feeder
- fully automatic test cycle
- special solid heating program

Datasheet:
12-1822.pdf

Ref.-No:
12-1822 (230V)
12-1823 (115V)



Flash-Point Testing with Closed Cup acc. ABEL (semi-automatic)

IP 170 - ISO 13 736 - NF T66-009
EN 924 - IP 491 - IP 492 - ISO 1516 - ISO 1523

Note! Abel-Pensky Model available upon request.

Datasheet:
12-1270.pdf

Ref.-No:
12-1270 (230 V)



Flash-Point Testing with Closed Cup acc. TAG (semi-automatic)

ASTM D 56 - FTM 791-1101 - JIS K2265

Datasheet:
12-1520.pdf

Ref.-No:
12-1520 (230V)
12-1521 (110V)



Flash-Point Testing with Open Cup acc. TAG (semi-automatic)

ASTM D 1310 - ASTM D 3143

Datasheet:
12-1540.pdf

Ref.-No:
12-1540 (230V)



Flash-Point Testing with Closed/Open Cup acc. RAPID Tester "Series 3" (automatic)

IP 523 - ASTM D 3828-B - ISO 3679 - JIS K2265
ASTM D 3278-B - UN Class 3 N - ASTM E 502

IP 524 - ASTM D 3828-A - ISO 3680
ASTM D 3278-A

ASTM D 4206 - ISO 9038 - BS 3900-A11 - UN Class 3 L.3 -
CFR 49-173.120

Datasheet:
12-1240.pdf

Ref.-No:
12-1240 (115/230V)
- closed cup
12-1241 (115/230V)
- open cup
12-1260 (115/230V)
- special



Flash-Point Testing with Closed Cup acc. PENSKEY-MARTENS - PMP 4 (semi-automatic)

ASTM D 93-A+B - JIS K2265 - FTM 791-1102 - FTM 141-4293 -
AASHTO T73 - AASHTO T172 - IP 34-A+B - ISO 2719-A+B
(EN 22 719 obs.)

ASTM D 6751

ISO 15 267 - BS 684-1.17

- 1-twist handling
- barometric pressure correction
- automatically controlled heating rate

Datasheet:
12-0595.pdf

Ref.-No:
12-0595 (115/230V)



Flash-Point Testing with Closed Cup acc. PENSKY-MARTENS - PM 4 (semi-automatic)

ASTM D 93-A+B - JIS K2265 - FTM 791-1102 - FTM 141-4293 -
AASHTO T73 - AASHTO T172 - IP 34-A+B - ISO 2719-A+B
(EN 22 719 obs.)

ASTM D 6751

ISO 15 267 - BS 684-1.17

- 1-twist handling
- gas & electric igniter
- manual adjusted heating rate

Datasheet:
12-1650.pdf

Ref.-No:
12-1650 (230V)



Flash-Point Testing with Closed Cup acc. PENSKY-MARTENS - PM (manual)

ASTM D 93-A - JIS K2265 - AASHTO T73-A - AASHTO T172-A -
FTM 791-1102 - FTM 141-4293 - IP 34-A - ISO 2719-A (EN 22 719 obs.)

ASTM D 6751

ISO 15 267 - BS 684-1.17

Optional IP 35

- adjustable gas burner (for all gas-types)
- hand operated stirrer
- open cup-option (IP 35)

Datasheet:
12-0601.pdf

Ref.-No:
12-0601 (Gas)



Flash-Point Testing with Closed Cup acc. PENSKY-MARTENS - Army Check (manual) (Field Test)

similar to international Pensky-Martens methods

Datasheet:
12-1000.pdf

Ref.-No:
12-1000 (230V)
12-1001 (115V)



Flash- & Fire-Point Testing with Open Cup acc. CLEVELAND - CL 5 (semi-automatic)

ASTM D 92 - JIS K2265 - AASHTO T48 -FTM 791-1103 - IP 36 -
ISO 2592

Superpave®

- re-cooling device
- electric heating with overheat protection
- gas ignition device

Datasheet:
12-2730.pdf

Ref.-No:
12-2730 (230V)
12-2731 (115V)



Flash- & Fire-Point Testing with Open Cup acc. CLEVELAND (manual)

ASTM D 92 - JIS K2265 - AASHTO T48 -FTM 791-1103 - IP 36 -
ISO 2592

Superpave®

- 3 interchangeable test flame jets
- adjustable gas burner (for all gas-types)
- gas ignition

Datasheet:
12-0701.pdf

Ref.-No:
12-0701 (Gas)



Auto-Ignition Temperature (semi-automatic)

ASTM E 659

Datasheet:
12-1802.pdf

Ref.-No:
12-1802 (230V)




Auto-Ignition Temperature of Solids acc. Grewer (semi-automatic)

VDI 2263 (in analogy to 79/831/EWG)

Datasheet:
12-1840.pdf

Ref.-No:
12-1840 (230V)





Water Vapor Content by Dew-Point (manual)

ASTM D 1142

Datasheet:
14-0004.pdf

Ref.-No:
14-0004
to 14-0006



Dryness by Freeze Valve (manual)

ASTM D 2713 - IP 395 - ISO 13 758

Datasheet:
14-0100.pdf

Ref.-No:
14-0100



Vapor Pressure (manual)

ASTM D 1267 - IP 161 - IP 410 - ISO 4256 - JIS K2240 - JIS K2258
- digital Manometer PA-REID available

Datasheet:
14-0015.pdf

Ref.-No:
14-0015



Corrosion by Copper Strip (manual)

ASTM D 1838 - ISO 2160 - ISO 6251 - IP 411 - JIS K2240

Datasheet:
14-0050.pdf

Ref.-No:
14-0050 (230V)
14-0047 (115V)




Volatility & Residue by Propane Residue (Mercury Freeze Method) (manual)

ASTM D 1837 - NF M41-012
ASTM D 2158 - IP 317

Datasheet:
14-0060.pdf

Ref.-No:
14-0060



Density or Relative Density by SCHILLING Effusimeter (manual)

IP 59-C (obs.)

Datasheet:
14-0090.pdf

Ref.-No:
14-0090



Density or Relative Density by Pressure Hydrometer (manual)

ASTM D 1657 - IP 235 - ISO 3993 (NF M41-008 obs.) - DIN 12 804 - DIN 51 757-5

Datasheet:
14-0110.pdf

Ref.-No:
14-0110 (230V)



Hydrogen Sulfide in Liquefied Petroleum (manual)

ASTM D 2420 - ISO 8819 (NF M41-011 obs.)

Datasheet:
14-0230.pdf

Ref.-No:
14-0230 (230V)
14-0231 (115V)



<p>Density Balance</p>	<p>see: Density</p>	<p>page 6</p>
<p>Flash-Point by TAG & Auto-Ignition</p>	<p>see: Flash-Point</p>	<p>page 12</p>
<p>GUM by Jet Evaporation (manual) ASTM D 381 - DIN 51 784 - FTM 791-3302 - ISO 6246 - IP 131 - IP 540 - JIS K2261 (BS 4348, EN 26 246 obs.) - highest level of user safety - all necessary auxiliary devices are available - complete re-calibration set is available</p>	<p>Datasheet: 13-0035.pdf Ref.-No: 13-0035 (230V) 13-0036 (115V)</p>	
<p>Aniline Point by Method A (Test Tube) (manual) ASTM D 611-A - IP 2-A - JIS K2256 - NF M07-021 ISO 2977-A - (DIN 51 775, DIN 51 787 obs.) - available with ASTM or DIN test tube</p>	<p>Datasheet: 13-0100.pdf Ref.-No: 13-0101 (ASTM) 13-0100 (DIN)</p>	
<p>Aniline Point by Method B (Thin Film) (manual) ASTM D 611-B - IP 2-B - JIS K2256 - NF M07-021 ISO 2977-B (DIN 51 775, DIN 51 787 obs.)</p>	<p>Datasheet: 13-0106.pdf Ref.-No: 13-0106 (230V) 13-0107 (115V)</p>	
<p>Aniline Point by Method D (U-Tube) (manual) ASTM D 611-D - IP 2-D - JIS K2256 - NF M07-021 ISO 2977-D (DIN 51 775, DIN 51 787 obs.)</p>	<p>Datasheet: 13-0102.pdf Ref.-No: 13-0102 (230V) 13-0103 (115V)</p>	
<p>Aniline Point by Method E (automatic) ASTM D 611-E - FTM 791-3601 - IP 2-E - JIS K2256 - NF M07-021 ISO 2977-E (DIN 51 775 - DIN 51 787 obs.) (Method E is the automatic procedure of methods A and B)</p>	<p>Datasheet: 13-0001.pdf Ref.-No: 13-0001 (230V)</p>	
<p>Fluorescent Indicator Adsorption - FIA (manual) ASTM D 1319 - FTM 791-3703 - IP 156 - ISO 3837 - JIS K2536 - NF M07-024 (EN 10 obs.) - 4-place - standard column test kit - precision bore column test kit</p>	<p>Datasheet: 13-0930.pdf Ref.-No: 13-0930 (230V) 13-0929 (115V)</p>	
<p>Depentanization by PONA (manual) ASTM D 2001</p>	<p>Datasheet: 13-1030.pdf Ref.-No: 13-1030 (230V) 13-1031 (115V)</p>	
<p>Vapor Pressure by REID (digital) ASTM D 323-A+B+D - ASTM D 4953-A - IP 69-A+B+D - ISO 3007-A+B+D - JIS K2258 ASTM D 323-C - IP 69-C - ISO 3007-C - JIS K2258 - digital manometer PA-Reid & baths available - liquid chamber (1 opening) for less than 180 kPa - liquid chamber (2 openings) for above than 180 kPa</p>	<p>Datasheet: 13-0201.pdf Ref.-No: 13-0201 (<180 bar) 13-0202 (>180 bar)</p>	



Vapor Pressure Analyzer - PetroVap II (automatic)

ASTM D 5191 - EN 13 016-1 - IP 394 - IP 481

Correlates with:

ASTM D 323 - JIS K2258 - ASTM D 4953 - ASTM D 5190 -
ASTM D 6378 - ISO 3007 - IP 69

Datasheet:
13-2620.pdf

Ref.-No:
13-2620 (115/230V)



Vapor Pressure Analyzer (automatic)

ASTM D 5191 - EN 13 016-1+2 - IP 394 - IP 409

Correlates with:

ASTM D 323 - JIS K2258 - ASTM D 2533 - ASTM D 4953 -
ASTM D 5190 - ASTM D 5482 - ASTM D 6378 - ISO 3007 - IP 69

Datasheet:
13-2600.pdf

Ref.-No:
13-2600 (105-250V)

page 21

Lead & Salt by Hot-Extraction - TEL

see: Constituents

page 21

Lead in Gasoline - Iodine Monochloride Method

see: Constituents



Freezing Point (manual)

ASTM D 2386 - IP 16 - ISO 3013 - JIS K2276 (DIN 51 421 obs.)

Datasheet:
13-0461.pdf

Ref.-No:
13-0461



Freezing Point (semi-automatic)

ASTM D 2386 - IP 16 - ISO 3013 - JIS K2276 (DIN 51 421 obs.)

- with 2 place cryostat

Datasheet:
13-1470.pdf

Ref.-No:
13-1470 (230V)
13-1471 (115V)



Freezing Point of Antifreeze Solutions (manual)

ASTM D 1177 - NF T78-102

Datasheet:
13-2460.pdf

Ref.-No:
13-2460 (230V)



Cold Filter Plugging Point - CFPP (manual)

ASTM D 6371 - EN 116 - IP 309 - JIS K2288
(DIN 51 428-1, SIS 155 122 obs.)

13-0420 CFPP - Dry-Ice Test Arrangement

Datasheet:
13-0420.pdf

Ref.-No:
13-0420



Cold Filter Plugging Point - CFPP I (manual)

ASTM D 6371 - EN 116 - EN 14214 - IP 309 - JIS K2288
(DIN 51 428-1, SIS 155 122 obs.)

13-2418 CFPP I - Cryostat with 1 Test-Chamber

Datasheet:
13-2418.pdf

Ref.-No:
13-2418 (230V)
13-2414 (115V)



Cold Filter Plugging Point - CFPP V (manual)

ASTM D 6371 - EN 116 - EN 14214 - IP 309 - JIS K2288
(DIN 51 428-1, SIS 155 122 obs.)

13-2419 CFPP V - Cryostat with 5 Test-Chambers

Datasheet:
13-2419.pdf

Ref.-No:
13-2419 (230V)

page 42

Refractive Index & Dispersion by Universal ABBE Refractometer by Modified ABBE Refractometer

see: Various

page 42

Refractive Index by High Performance Refractometer - ABBEMAT by Hand Refractometer - Eclipse by Refractometer - RFM 340

see: Various



Conductivity by Portable Method - MLA
(manual) (Field Test)

ASTM D 2624 - DIN 51 412-2 - ISO 6297 - IP 274 - JIS K2276

Datasheet:
13-3180.pdf

Ref.-No:
13-3180 (Battery)



Conductivity by Portable Method - DCM
(manual) (Field Test)

ASTM D 2624 - DIN 51 412-2 - ISO 6297 - IP 274 - JIS K2276

Datasheet:
13-0195.pdf

Ref.-No:
13-0195 (Battery)



Conductivity by Lab Method - Precision (automatic)

ASTM D 4308

Datasheet:
13-1410.pdf

Ref.-No:
13-1410 (Battery)

Water Separation Characteristics
by Micro-Separometer - MSEP (Field Test)

ASTM D 3948 - ASTM D 4860 - JIS K2276
ASTM D 7724 - DIN 51403-2

- with petrotest's exclusive "Safe-T-Close"-System

Datasheet:
13-0391.pdf

Ref.-No:
13-0391



Water Reaction by Mixing Cylinder (manual)

ASTM D 1094

Datasheet:
13-0401.pdf

Ref.-No:
13-0401



Oxidation Stability (RSSOT) - PetroOXY (automatic)

ASTM D 7525 (Gasoline)
ASTM D 7545 - EN 16091 (Diesel, Biodiesel (FAME) & Blends)
Correlates with:
Gasoline Method: ASTM D 525 - IP 40 - ISO 7536
Diesel Method: ASTM D 2274 - IP 388 ISO 12205; ASTM D 5304
FAME Method: EN 14112 - EN 15751
Grease Method: ASTM D 942 - DIN 51808 - IP 142
Oil Method: ASTM D 943 - ISO 4263; ASTM D 2272 - IP 229

Datasheet:
13-3002.pdf
13-3006.pdf

Ref.-No:
13-3002 (115/230V, EU-plug)
13-3003 (115/230V, US-plug)



PetroOXY "Block-Cooling"

- manual re-cooling between tests by chilled cooling block

PetroOXY "Active Cool"

- automatic re-cooling between tests by peltier elements

Ref.-No:
13-3006 (115/230V, EU-plug)
13-3007 (115/230V, US-plug)



Oxidation Stability - OBA (semi-automatic)

ASTM D 525 - FTM 791-3352 - IP 40 - ISO 7536 - JIS K2287
(BS 4347, DIN 51 780 obs.)

ASTM D 873 - FTM 791-3354 - IP 138 - JIS K2276
(BS 4456, DIN 51 799 obs.)

- 4-place Liquid or Dry-Heat baths with vessel-holding option
- safety approved oxidation vessels available
- digital Manometer

Datasheet:
13-1660.pdf

Ref.-No:
13-1660 (230V)
13-1661 (115V)



Oxidation Stability by
Liquid Bath Test Arrangement - OBA 1

Oxidation Stability by
Dry-Heat Bath Test Arrangement - OBA 1/T

Ref.-No:
13-1662 (230V)
13-1663 (115V)





Oxidation Stability by Liquid Bath Safety Test Arrangement - OBA 2

Datasheet:
13-1664.pdf

Ref.-No:
13-1664 (230V)
13-1665 (115V)

Oxidation Stability by Dry-Heat Bath Safety Test Arrangement - OBA 2/T

Ref.-No:
13-1666 (230V)
13-1667 (115V)



Smoke Point by Lamp Method (manual)

ASTM D 1322 - FTM 791-2107 - IP 57 - ISO 3014 - JIS K2537 -
NF M07-028 (DIN 51 406 obs.)

Datasheet:
13-0691.pdf

Ref.-No:
13-0691

page 23

Color by SAYBOLT-Scale Chromometer

see: Oil



Corrosion by Copper Strip Tarnish Test (manual)

ASTM D 130 - ASTM D 4048 - DIN 51 811 - FTM 791-5325 - IP 112 -
IP 154 - ISO 2160 - JIS K2220 - JIS K2513
(BS 4351, DIN 51 759-1 obs.)

Datasheet:
13-0700.pdf

Ref.-No:
13-0700 (230V)
13-0699 (115V)



Corrosion by Silver Strip Test (manual)

IP 227 (obs.)

Datasheet:
13-1000.pdf

Ref.-No:
13-1000 (230V)



Corrosion by Copper Strip Test (manual)

ASTM D 849

13-1020 with 1 test-place

13-1028 with 4 test-places

Datasheet:
13-1020.pdf

Ref.-No:
13-1020 (230V)
13-1019 (115V)

Ref.-No:
13-1028 (230V)
13-1029 (115V)



Coking Tendency by CONRADSON Carbon Residue (manual)

ASTM D 189 - DIN 51 551-1 - FTM 791-5001 - IP 13 - ISO 6615 -
JIS K2270 (BS 4380, NF T60-116 obs.)

13-0841 for MUNICIPAL gas supply
13-0842 for BUTANE/PROPANE gas supply
13-0843 for NATURAL gas supply

Datasheet:
13-0841.pdf

Ref.-No:
13-0841 (Gas)
13-0842 (Gas)
13-0843 (Gas)



Coking Tendency by Micro Carbon Residue - MCR-210 (manual)

ASTM D 4530 - IP 398 - ISO 10 370 - JIS K2270

Datasheet:
13-2900.pdf

Ref.-No:
13-2900 (230V)






Coking Tendency by RAMSBOTTOM Carbon Residue - RCRT (manual)

ASTM D 524 - FTM 791-5002 - IP 14 - ISO 4262 - NF T60-117
(BS 4451 obs.)

Datasheet:
13-0870.pdf

Ref.-No:
13-0870 (230V)
13-0871 (115V)



Particulates & Filter Membrane Color Ratings in Aviation Fuel by In-Line Sampler - AFCON	see: Constituents	page 21
Particulates & Free Water	see: Constituents	page 21
Particulate Contamination in Aviation Fuel by Laboratory Filtration ASTM D 5452 - IP 423	Datasheet: 16-1960.pdf Ref.-No: 16-1960	
Heat of Combustion by Calorimeter - C 2000 (automatic) ASTM D 240 - ASTM D 4809 - ASTM D 5468 - DIN 51 900-1-3 - ISO 1928 - JIS K2279 (ASTM D 1989 obs.)	Datasheet: 13-1880.pdf Ref.-No: 13-1880 (230V)	
Water & Sediment by Centrifuge	see: Constituents	page 22
Sediment by Extraction	see: Constituents	page 22
Total Sediment	see: Constituents	page 22
Salts by Conductivity (electrometric) - Salinometer	see: Constituents	page 21
Asphaltenes by Extraction	see: Constituents	page 21
Mercaptan by Doctor-Test	see: Constituents	page 21
Stability - Thermal for Jet Fuel (manual) ASTM D 3241 - FTM 791-2506 - IP 323 - ISO 6249 (ASTM D 1661 obs.)	Datasheet: 13-0272.pdf Ref.-No: 13-0272	
Sulfur & Chlorine by WICKBOLD Combustion	see: Constituents	page 20



Sulfur & Chlorine by WICKBOLD Combustion (manual)

ASTM D 2384 - ASTM D 2784 - DIN 51 408-1 - EN 24 260 -
FD M07-008 - IP 243 - ISO 4260 - JIS K2240 - JIS K2541-1

Datasheet:
13-0670.pdf

Ref.-No:
13-0670 (230V)
13-0671 (115V)



Sulfur & Chlorine by GROTE-KREKELER (manual)

DIN 51 400-2 - JIS K2541

Datasheet:
13-0756.pdf

Ref.-No:
13-0756 (230V)



Sulfur & Chlorine by High Temperature (manual)

ASTM D 1552 - DIN 51 724-1 - DIN 51 727

Datasheet:
13-1160.pdf

Ref.-No:
13-1160 (230V)
13-1161 (115V)




Sulfur by SCHÖNIGER (Thorin-Titration) (manual)

DIN 51 400-3 (similar IP 242)

Datasheet:
13-1720.pdf

Ref.-No:
13-1720



Sulfur by Raney-Nickel (manual)

DIN 51 400-8 (obs.)

Datasheet:
13-0475.pdf

Ref.-No:
13-0475 (230V)



Sulfur (Active) (manual)

ASTM D 1662

Datasheet:
16-0084.pdf

Ref.-No:
16-0084 (230V)
16-0085 (115V)



Sulfonation Number by Shaking Machine (manual)

FTM 791-1040

Datasheet:
13-0803.pdf

Ref.-No:
13-0803 (230V)



Unsulfonated Residue of Plant-Spray-Oil by Combined Shaking Machine (manual)

ASTM D 483

Datasheet:
13-0920.pdf

Ref.-No:
13-0920 (230V)



Mercaptan Sulfur by pH-Meter (manual)

ASTM D 3227 - IP 342 - ISO 3012 - JIS K2276 (NF M07-022 obs.)

Datasheet:
13-1560.pdf

Ref.-No:
13-1560 (230V)

Mercaptan by Doctor-Test (manual)

ASTM D 4952 - DIN 12 685-2 - IP 30 - ISO 5275 - JIS K2276

Datasheet:
16-0431.pdf

Ref.-No:
16-0431



Asphaltenes by Extraction (manual)

ASTM D 6560 - DIN 51 595 - IP 143

Datasheet:
13-0170.pdf

Ref.-No:
13-0170 (230V)
13-0169 (115V)



Fluorescent Indicator Adsorption - FIA

see: Fuel

page 15

Depentanization by PONA

see: Fuel

page 15

**Salts by Conductivity (electrometric)
- Salinometer SCO 1 (manual)**

ASTM D 3230 - IP 265 - JIS K2601

- laboratory bench model

Datasheet:
13-0497.pdf

Ref.-No:
13-0497 (230V)



**Salts by Conductivity (electrometric)
- Salinometer Portable (manual) (Field Test)**

ASTM D 3230 - IP 265 - JIS K2601

Datasheet:
13-1490.pdf

Ref.-No:
13-1490 (230V)
13-1491 (115V)



Lead & Salt by Hot-Extraction - TEL (manual)

IP 77 - DIN 51 576

FTM 791-5501 - IP 182 (ASTM D 526, ASTM D 2547, DIN 51 769-1, IP 96, IP 248, ISO 2083 obs.)

Datasheet:
13-0305.pdf

Ref.-No:
13-0305 (220V)



Lead in Gasoline - Iodine Monochloride Method (manual)

ASTM D 3341

Datasheet:
13-2340.pdf

Ref.-No:
13-2340



Particulates & Filter Membrane Color Ratings in Aviation Fuel by In-Line Sampler - AFCON (manual)

ASTM D 2276 - IP 216 - FTM 791-3008 - JIS K2276

Datasheet:
13-1980.pdf

Ref.-No:
13-1980



Particulates & Free Water (manual)

ASTM D 4176

Datasheet:
13-2350.pdf

Ref.-No:
13-2350



Particulates by Diaphragm Filtration (manual)

EN 12 662 - IP 440

Datasheet:
16-0970.pdf

Ref.-No:
16-0970



**Sediment, Total (manual)**

ASTM D 4870 - IP 375 - IP 390-A - ISO 10 307-1+2

Datasheet:
13-2300.pdfRef.-No:
13-2300 (230V)
13-2301 (115V)**Sediment by Extraction (manual)**

ASTM D 473 - FTM 791-3002 - IP 53 - ISO 3735 (BS 4382 obs.)

Datasheet:
13-0321.pdfRef.-No:
13-0321**Water & Sediment by Centrifuge - PetroCEN
(semi-automatic)**Lubricating Oil Tests:

ASTM D 91 - ASTM D 2273 - JIS K2503

Datasheet:
13-2140_ASTM D 91.pdf
13-2140_ASTM D 2273.pdfRef.-No:
13-2140 (230V)
13-2141 (120V)**Water & Sediment by Centrifuge - PetroCEN
(semi-automatic)**Lubricating Oil Tests:

ASTM D 893

Amount of Water and Sediment:

API 2542 - AASHTO T55 - JIS K2601 (ASTM D 96, BS 4385 obs.)

ASTM D 4007 (IP 359 obs.)

ASTM D 1796 - API 2548 (DIN 51 793 - IP 75 obs.)

ASTM D 2709

ISO 3734

ISO 9030

Datasheet:
13-2160_ASTM D 96.pdf
13-2160_ASTM D 893.pdf
13-2160_ASTM D 1796.pdf
13-2160_ASTM D 2709.pdf
13-2160_ASTM D 2711.pdf
13-2160_ASTM D 4007.pdf
13-2160_ISO 3734.pdf
13-2160_ISO 9030.pdfRef.-No:
13-2160 (230V)
13-2161 (120V)
- heatable < +70°C

**Demulsification Number (manual)**

IP 19 (ASTM D 1935 obs.)

Datasheet:
16-0140.pdfRef.-No:
16-0140 (230V)
16-0141 (115V)**Demulsibility Characteristics (manual)**

ASTM D 2711

Datasheet:
16-1230.pdfRef.-No:
16-1230 (230V)**Demulsibility Characteristics by HERSCHEL Stirring Method (semi-automatic)**ASTM D 1401 - IP 412 - ISO 6614 - JIS K2207 - JIS K2520
(DIN 51 599, FTM 791-380 obs.)

- 7 place full-view bath
- variable stirrer speed with digital display
- automatic timer stop and audible alarm

Datasheet:
16-2100.pdfRef.-No:
16-2100 (230V)
16-2101 (115V)**Color by SAYBOLT-Scale Chromometer (manual)**

ASTM D 156 - DIN 51 411 - FTM 791-101 - JIS K2580 - NF M07-003

with scales for: Saybolt Color

Datasheet:
16-0180.pdfRef.-No:
16-0180 (230V)
16-0179 (115V)**Color by ASTM-Scale Comparator (manual)**ASTM D 1500 - IP 196 - ISO 2049 - JIS K2580
(DIN 51 578, FTM 791-102, NF T60-104 obs.)

with scales for: ASTM Color

Datasheet:
16-0210.pdfRef.-No:
16-0210 (115/230V)**Color by Colorimeter - PFX 195/2 (automatic)**with scales for: Saybolt Color, ASTM Color, CIE Values, Spectral Data,
Pt-Co / Hazen / APHAASTM D 6045
ASTM D 156 - DIN 51 411 - JIS K2580 - NF M07-003
ASTM D 1500 - IP 196 - ISO 2049 - JIS K2580
(DIN 51 578, FTM 791-102, NF T60-104 obs.)
ASTM D 1209 - ASTM D 6166
(ASTM E 308 obs.)optional scales for: Gardner

ASTM D 1544

Datasheet:
16-2194.pdf
Ref.-No:
16-2194 (115/230V)**Color by Colorimeter - PFX 880/P (automatic)**with scales for: Tristimulus Method (Automatic Saybolt Color &
ASTM Color Scale), CIE Values, Spectral DataASTM D 6045
ASTM D 156 - DIN 51 411 - JIS K2580 - NF M07-003
ASTM D 1500 - IP 196 - ISO 2049 - JIS K2580
(DIN 51 578, FTM 791-102, NF T60-104 obs.)
(ASTM E 308 obs.)optional scales for: Pt-Co / Hazen / APHA, GardnerASTM D 1209 - ASTM D 6166
ASTM D 1544Datasheet:
16-2175.pdf
Ref.-No:
16-2175 (115/230V)
16-2176 (115/230V)

16-2175 Colorimeter PFX 880/P

16-2176 Colorimeter PFX 880/P-H with heated sample chamber



Color by Colorimeter - PFX 995/P (automatic)



with scales for: Tristimulus Method (Automatic Saybolt Color & ASTM Color Scale), CIE Values, Spectral Data, Pt-Co / Hazen / APHA, Gardner, Lovibond RYBN, IP Units

ASTM D 6045
ASTM D 156 - DIN 51 411 - JIS K2580 - NF M 07-003
ASTM D 1500 - IP 196 - ISO 2049
(DIN 51 578, FTM 791-102, NF T 60-104 obs.)
ASTM D 1209 - ASTM D 6166
ASTM D 1544
IP 17-A+B
(ASTM E 308 obs.)

16-2190 Colorimeter PFX 995/P
16-2191 Colorimeter PFX 995/P-H with heated sample chamber

Datasheet:
16-2190.pdf

Ref.-No:
16-2190 (115/230V)
16-2191 (115/230V)



Tensiometer - Processor (automatic)

ASTM D 971 - ASTM D 1331 - ASTM D 1417 - EN 14 370 - ISO 1409 - ISO 6295 - ISO 6889 - JIS K2241 (ASTM D 1590, DIN 53 914 obs.)

Lecomte du Noüy Ring & Wilhelmy Plate Method

Datasheet:
16-1750.pdf

Ref.-No:
16-1750 (115/230V)



Tensiometer - Scale (manual)

ASTM D 971 - EN 14 370 - JIS K2241 (DIN 53 914 obs.)

Lecomte du Noüy Ring Method

Datasheet:
16-1755.pdf

Ref.-No:
16-1755



Total Dirt by Refrigerated Centrifuge (semi-automatic)

DIN 51 365

Datasheet:
16-1870.pdf

Ref.-No:
16-1870 (230V)



Fluidity by U-Tube Method (manual)

DIN 51 568

Datasheet:
16-0252.pdf

Ref.-No:
16-0252 (230V)
16-0253 (115V)



Flock Point (manual)

DIN 51 351

Datasheet:
16-0390.pdf

Ref.-No:
16-0390 (230V)
16-0389 (115V)



Resistance by Philipp Test (manual)

DIN 51 593 (obs.)

Datasheet:
16-0401.pdf

Ref.-No:
16-0401 (230V)
16-0402 (115V)



Solidity & Turbidity Point (manual)

DIN 51 583 (obs.)

Datasheet:
16-0301.pdf

Ref.-No:
16-0301

**Cloud and Pour Point - CAPP 1 & CAPP 4 (manual)**

ASTM D 97 - ASTM D 2500 - ASTM D 5853 - EN 23 015 -
FTM 791-201 - IP 15 - IP 219 - IP 441 - ISO 3015 - ISO 3016 -
NF T60-105 - JIS K2269 - JIS K2601 - (BS 4458, DIN 51 597 obs.)

- 1-place type (16-0352)
- 4-place type (16-2348 & 16-2349)
- both types with circulation coil for external cooling

Datasheet:
16-0352.pdf

Ref.-No:
16-0352
16-2348 (230V)
16-2349 (115V)



CAPP 1

CAPP 4

Cloud and Pour Point - CAPP I (semi-automatic)

ASTM D 97 - ASTM D 2500 - ASTM D 5853 - EN 23 015 -
FTM 791-201 - IP 15 - IP 219 - IP 441 - ISO 3015 - ISO 3016 -
NF T60-105 - JIS K2269 - JIS K2601 - (BS 4458, DIN 51 597 obs.)

1 chamber with 4 test-positions & metal-block cryostat

Datasheet:
16-2372.pdf

Ref.-No:
16-2372 (230V)
16-2373 (115V)

**Cloud and Pour Point - CAPP V (semi-automatic)**

ASTM D 97 - ASTM D 2500 - ASTM D 5853 - EN 23 015 -
FTM 791-201 - IP 15 - IP 219 - IP 441 - ISO 3015 - ISO 3016 -
NF T60-105 - JIS K2269 - JIS K2601 - (BS 4458, DIN 51 597 obs.)

5 chambers with 4 test-positions each & metal-block cryostat

Datasheet:
16-2360.pdf

Ref.-No:
16-2360 (230V)
16-2363 (115V)

**Shear Stability by Diesel Injector Rig (manual)**

ASTM D 2603 - ASTM D 6278 - ASTM D 7109 - CEC L-14-A-95 -
DIN 51 382 - IP 294

Datasheet:
16-0320.pdf

Ref.-No:
16-0320 (380V)
16-1383 (3x 220V)

**Shear Stability of Polymer-Containing Oils and Hydraulic Fluids - Sonic (semi-automatic)**

ASTM D 2603 - ASTM D 5621

Datasheet:
16-2410.pdf

Ref.-No:
16-2410 (230V)

**Viscosity, High Temperature & High Shear Rate by Tapered-Plug Viscometer - HTHS (semi-automatic)**

ASTM D 4741 - CEC L-36-A-90

Datasheet:
16-1660.pdf

Ref.-No:
16-1660 (230V)

**Viscosity at Low Temperature by Brookfield (semi-automatic)**

ASTM D 2983
DIN 51 398 - IP 267-B

Datasheet:
16-1152.pdf

Ref.-No:
16-1152 (230V)

**Viscosity at Low Temperature by Brookfield (automatic)**

ASTM D 2983

available Viscometer: LVDV-I+ or LVDL-II+ Pro

with Bath Temperature Range: +30 to -40 °C
or -----

with Bath Temperature Range: +30 to -70 °C

Datasheet:
16-2500.pdf

Ref.-No:
16-2500 (220V, 50 Hz)
16-2501 (220V, 60 Hz)
16-2502 (110V, 60 Hz)

16-2503 (220V, 50 Hz)
16-2504 (220V, 60 Hz)





Viscosity at Low Temperature (Borderline Pumping Temp.) by Mini-Rotary Viscometer - MPV/TP-1 (automatic)

ASTM D 4684 - ASTM D 6896 - ASTM D 3829

by Constant Shear Stress Viscometer (automatic)

ASTM D 6821

Datasheet:
16-2430.pdf

Ref.-No:
16-2430 (220V, 50 Hz)
16-2431 (220V, 60 Hz)
16-2432 (110V, 60 Hz)



Oil Content of Water Mixed Cooling Lubricants by Acidification & Saponification (manual)

DIN 51 367 - DIN 51 385
DIN 51 368 - IP 137

Datasheet:
16-0070.pdf

Ref.-No:
16-0070 (230V)
16-0072 (230V)



Stability by Oil / Layer Separation (manual)

DIN 51 346
IP 263 (similar)

16-0074 with emulsion stirring unit
16-0082 with magnetic stirrer unit

Datasheet:
16-0074.pdf

Ref.-No:
16-0074 (230V)
16-0082 (230V)



Corrosion Preventing Characteristics by HERBERT Corrosion Test (manual)

DIN 51 360-1 - IP 125

Datasheet:
16-0281.pdf

Ref.-No:
16-0281 (230V)
16-0282 (115V)



Corrosion Preventing Characteristics by Chip & Filterpaper Test (manual)

DIN 51 360-2

Datasheet:
16-0286.pdf

Ref.-No:
16-0286 (230V)
16-0287 (115V)

page 22

Insoluble Matters & Total Contamination by Diaphragm Filtration

see: Constituents

page 20

Sulfur (Active)

see: Constituents



pH-Value by pH-Meter (manual)

DIN 51 369 (partially for ASTM D 1121, ASTM D 1287)

Datasheet:
16-0295.pdf

Ref.-No:
16-0295



Ash by Muffle Furnace (manual)

ASTM D 482 - ASTM D 874 - ASTM D 1026 - ASTM D 4422 -
DIN 51 352-2 - DIN 51 575 - FTM 791-5421 - FTM 791-5422 -
IP 4 - IP 163 - ISO 3987 - ISO 6245 - JIS K2272

Datasheet:
16-0001.pdf

Ref.-No:
16-0001 (230V)



Aging Stability by Baader (manual)

DIN 51 554-1+2

Datasheet:
16-0750.pdf

Ref.-No:
16-0750 (230V)
16-0749 (115V)



Gasoline Diluent of Engine Oil by Distillation (manual)

ASTM D 322 - IP 23

16-0683 flask heating on electric heater
16-0683 flask heating with heating mantle

Datasheet:
16-0683.pdf

Ref.-No:
16-0683 (230V)
16-1317 (230V)



Oxidation Characteristics by TOST/TOO Dry-Heat - ASTM (manual)

4-place dry-heat bath: for tests up to 230°C

ASTM D 943 - ISO 4263 - ISO 12 205 - JIS K2514 (BS 4388 obs.)

Datasheet:
16-0586_ASTM D 943.pdf
16-0586_ASTM D 4310.pdf
16-0586_DIN 51587.pdf
16-0586_ISO 12205.pdf

Ref.-No:
16-1416 (230V)
16-1417 (115V)



Oxidation Characteristics by TOST/TOO Dry-Heat - DIN (manual)

4-place dry-heat bath: for DIN-tests up to 230°C

(DIN 51 394 obs.) - FTM 791-5308

Datasheet:
16-0545.pdf

Ref.-No:
16-1418 (230V)
16-1419 (115V)

Oxidation Characteristics by TOST/TOO Liquid - ASTM D 943 (manual)

4-place liquid bath: for tests up to 175°C

ASTM D 943 - ASTM D 2274 - ASTM D 4310 - IP 157 - IP 388 -
ISO 12 205 - DIN 51 587 - (DIN 51 394 obs.) - FTM 791-5308 -
JIS C2101 - JIS K2514 (BS 4388 obs.)

8-place liquid bath: for tests up to 175°C

ASTM D 943 - ASTM D 2274 - ASTM D 4310 - IP 157 - IP 388 -
ISO 12 205 - DIN 51 587 - (DIN 51 394 obs.) - FTM 791-5308 -
JIS C2101 - JIS K2514 (BS 4388 obs.)

Datasheet:
16-0586_ASTM D 943.pdf
16-0586_ASTM D 4310.pdf
16-0586_DIN 51587.pdf
16-0586_ISO 12205.pdf
16-0545.pdf
16-0586_ASTM D 2274.pdf

Ref.-No:
16-0583 (230V)
16-0584 (115V)

Ref.-No:
16-0581 (230V)
16-0458 (115V)



Oxidation Characteristics by TOST/TOO Liquid - ASTM D 2893 (manual)

4-place liquid bath: for tests up to 175°C

ASTM D 2893

Datasheet:
16-0539.pdf

Ref.-No:
16-1420 (230V)
16-1421 (115V)

Oxidation Stability (RSSOT) - PetroOXY (automatic)

see: Fuel

page 17

Oxidation Stability & Corrosiveness (manual)

ASTM D 4636

- aluminum block heater

Datasheet:
25-0716.pdf

Ref.-No:
25-0716 (230V)



Oxidation Stability by BAM 200 Liquid (manual)

BAM 200 Liquid: flowmeter range 0.08...2 l/h

ASTM D 2440 - DIN 51 373 - DIN 51 538 - EN 61 125-A,B,C -
IP 280 - IP 306 - IP 307 - IP 335 - ISO 7624 - JIS C2101

- 8-place liquid bath with low level and temperature excess protection

Datasheet:
16-1500_ASTM D 2440.pdf
16-1500_DIN 51373.pdf
16-1500_DIN 51538.pdf
16-1500_EN 61 125.pdf
16-1500_IP 280.pdf
16-1500_IP 306.pdf

Ref.-No:
16-1500 (230V)
16-1501 (115V)



Oxidation Stability by BAM 250 Liquid (manual)

BAM 250 Liquid: flowmeter range 1.6...16 l/h

DIN 51 352-1+2 - IP 48 - ISO 6617-1+2

- 8-place liquid bath with low level and temperature excess protection

Datasheet:
16-1502_DIN 51352-1.pdf
16-1502_DIN 51352-2.pdf
16-1502_IP 48.pdf

Ref.-No:
16-1502 (230V)
16-1503 (115V)



Oxidation Stability by BAM 200 Dry-Heat (manual)

BAM-200 Dry: flowmeter range 0.08...2 l/h (Tube Ø 28)
DIN 51 373

BAM-200 Dry: flowmeter range 0.08...2 l/h (Tube Ø 26)
DIN 51 373

- 8-place dry-heat bath with aluminum block heater
- Option: available with 12-test places

Datasheet:
16-1500_DIN 51373.pdf

Ref.-No:
16-1519 (230V)
16-1520 (115V)

Datasheet:
16-1500_ASTM D 2440.pdf
16-1500_DIN 51538.pdf
16-1500_EN 61 125.pdf
16-1500_IP 280.pdf
16-1500_IP 306.pdf

Ref.-No:
16-1566 (230V)
16-1568 (115V)

Oxidation Stability by BAM 250 Dry-Heat (manual)

BAM-250 Dry: flowmeter range 1.6...16 l/h (Tube Ø 28)
DIN 51 352-1+2 - IP 48 - ISO 6617-1+2

BAM-250 Dry: flowmeter range 1.6...16 l/h (Tube Ø 26)
DIN 51 352-1+2 - IP 48 - ISO 6617-1+2

- 8-place dry-heat bath with aluminum block heater
- Option: available with 12-test places

Datasheet:
16-1502_DIN 51352-1.pdf
16-1502_DIN 51352-2.pdf

Ref.-No:
16-1522 (230V)
16-1523 (115V)

Datasheet:
16-1502_IP 48.pdf

Ref.-No:
16-1569 (230V)
16-1570 (115V)



Oxidation Stability by Rotating Bomb - RBOT (automatic)

ASTM D 2112 - ASTM D 2272 - IP 229 - JIS K2514 - ASTM D 4742

Datasheet:
16-1980.pdf

Ref.-No:
16-1980 (230V)



Oxidation Stability by Rotating Bomb - RBOT (semi-automatic)

ASTM D 2112 - ASTM D 2272 - IP 229 - JIS K2514 - ASTM D 4742

- 2-place bath
- digital manometers PA4-RBOT
- pressure indication in psi, kPa, bar, mWc

Datasheet:
16-1990.pdf

Ref.-No:
16-1990 (230V)
16-1991 (115V)



Rust Preventing Characteristics - TOR (semi-automatic)

ASTM D 665 - BS 7120 - FTM 791-4011 - IP 135 - ISO 7120 - JIS K2510 - NAC TM0172-2001

ASTM D 3603 - ASTM D 5534 - FTM 791-5315

Datasheet:
16-0724.pdf

Ref.-No:
16-0724 (230V)
16-0725 (115V)



Rust Preventing by Chip & Filter (manual)

DIN 51 360-2
DIN 51 360-1 - IP 125 - IP 287

Datasheet:
16-0249.pdf

Ref.-No:
16-0249
16-0242



Rust Protection by Humidity Cabinet (manual)

ASTM D 1748 - JIS K2246

Datasheet:
16-3250.pdf

Ref.-No:
16-3250 (230V)



Compatibility of Elastomer or Rubber with Lubricating Fluids & Greases

ASTM D 4289

Datasheet:
17-0780.pdf

Ref.-No:
17-0780 (230V)

Air Release Value by Impinger (semi-automatic)

ASTM D 3427 - IP 313 - ISO 9120

Datasheet:
16-1200.pdf

Ref.-No:
16-1200 (230V)
16-1201 (115V)



Foaming Characteristics - Quadruple (manual)

ASTM D 892 (sequence I-III) - IP 146 - ISO 6247 - JIS K2518

Optional:

ASTM D 6082 (sequence IV)

Twin Foam-Tester

- 4-places in each bath
- full-view baths with digital-thermostats
- each bath to be maintained individually at 24, 49, 93.5 or 150 °C

Single Foam-Tester

- 4-place
- full-view bath with digital-thermostat
- to be maintained at 24, 49, 93.5 or 150 °C

Datasheet:
16-1470.pdf

Ref.-No:
16-1470 (230V)
16-1471 (115V)



Foaming Tendencies of Engine Coolants (manual)

ASTM D 1881

Datasheet:
16-0690.pdf

Ref.-No:
16-0690 (230V)
16-0689 (115V)



Corrosion of Engine Coolants (manual)

ASTM D 1384 - JIS K2234

Datasheet:
16-0955.pdf

Ref.-No:
16-0955 (230V)



Evaporation Loss by NOACK - VP 4000 (automatic)

Not yet approved by Standards Organisations but in analogy to:
ASTM D 5800-B - CEC L-40-A-93 - IP 421-B

wood metal free!

Datasheet:
16-2010.pdf

Ref.-No:
16-2010 (230V)
16-2011 (115V)



Evaporation Loss of Lubricants (manual)

ASTM D 972 - ASTM D 2878 - FTM 791-351 - IP 183 - JIS K2220

Datasheet:
16-0771.pdf

Ref.-No:
16-0771 (230V)
16-0772 (115V)



Hydrolytic Stability by Beverage Bottle (manual)

ASTM D 2619 - DIN 51 348

Datasheet:
99-0117.pdf

Ref.-No:
99-0117 (230V)



Dielectric Breakdown Voltage (automatic)

IEC 156 Fig.I - EN 60156 - IP 295 - JIS C 2101 - IEC 60156
ASTM D 1816 - IEC 156 Fig.II - VDE 0370-5 - EN 60156 - IP 295
ASTM D 877

Voltage ranges: 0 ... 60 kVrms or 0 ... 80 kVrms or 0 ... 100 kVrms
Electrode shapes: spherical / ball or mushroom / VDE / dome or disk

Datasheet:
16-2740.pdf

Ref.-No:
16-2740 (60 kV)
16-2745 (80 kV)
16-2750 (100 kV)





Neutralization Number by Color-Indicator Titration (Acid & Base No.) - TAN & TBN (manual)

ASTM D 974 - DIN 51 558-1+2 - FTM 791-5105 - IP 139 - JIS C2101 - JIS K2501 (BS 2634 obs.)

Datasheet: 17-0600.pdf

Ref.-No: **17-0600**



Neutralization Number by Potentiometric Titration (Acid & Total Base No.) - TAN & TBN (manual)

ASTM D 664 - ASTM D 4739 (alternative) - FTM 791-5106 - IP 177 - JIS K2501 (BS 4457, DIN 51 809-2 obs.)
ASTM D 2896 - IP 276 - ISO 3771 - JIS K2501

Datasheet: 17-0695.pdf

Ref.-No: **17-0695** (230V)
17-0696 (115V)



Acid Number by Potentiometric Titration - TAN (automatic)

ASTM D 664 - IP 177

Datasheet: 17-1800.pdf

Ref.-No: **17-1800** (110-240V)

page 42

Refractive Index by Refractometer

see: Fuel



Coking Tendency by Panel Coking (manual)

FTM 791-3462

Datasheet: 99-0300.pdf

Ref.-No: **99-0300** (220/240V)

page 18

Coking Tendency by Carbon Residue

see: Fuel

page 41

Water Content

see: Various



Lube Oil - VARIO Check (Field Test)

testing engine fluids on:
Viscosity, Water in Oil, Base No, Salt Water Contamination, Insolubles

Datasheet: 40-0160.pdf

Ref.-No: **40-0160**



Viscosity by Comparator - Viscosity Check (Field Test)

Datasheet: 40-0150.pdf

Ref.-No: **40-0150**



Viscosity by Falling Ball - JUNG Check (Field Test)

Datasheet: 40-0140.pdf

Ref.-No: **40-0140**



Water in Oil - WIO Check (Field Test)

Datasheet: 40-0120.pdf

Ref.-No: **40-0120**



PCB-Content - PCB Check (Field Test)

USEPA SW-846-9079
(CLOR-N-OIL)

Datasheet: 13-1650.pdf

Ref.-No: **13-1650**


Friction & Wear by REICHERT M2 (manual)

Datasheet:
15-0035.pdf
Ref.-No:
15-0035 (230V)


Lubrication & Load Capability by BRUGGER (manual)

DIN 51 347

Datasheet:
15-0070.pdf
Ref.-No:
15-0070 (380V)


Wear Characteristics (Extreme Pressure & Wear Properties, Shear Stability) by SHELL - Four-Ball Machine (manual)

ASTM D 2266 - FTM 791-6514
ASTM D 2596
ASTM D 2783 - JIS K2519
ASTM D 4172
ASTM D 5183
IP 239
DIN 51 350-1-6
FTM 791-6503
CEC L-45-T-99
ISO 20 623

With additional equipment:

ASTM D 2266 - ASTM D 2783 - ASTM D 4172 - CEC L-45-T-99 -
DIN 51 350-6 - PV 1417 - PV 1437 - PV 1444

Datasheet:
15-1000.pdf
Ref.-No:
15-1000 (380V)


Wear Characteristics (Extreme Pressure & Wear Properties, Load Carrying Capacity) by FALEX Pin & Vee Block (manual)

ASTM D 2625
ASTM D 2670
ASTM D 3233

Datasheet:
15-0450.pdf
Ref.-No:
15-0450 (230V)


Wear Characteristics (Extreme Pressure Properties, Load-Carrying Capacity) by TIMKEN-Test (manual)

ASTM D 2509 - IP 326 - JIS K2220
ASTM D 2782 - IP 240 - JIS K2519

Datasheet:
15-0602.pdf
Ref.-No:
15-0602 (230V)
15-0603 (115V)


Wear Characteristics (Scuffing Load Capacity & Shear Stability) by FZG Gear Oil Rig (manual)

ASTM D 4998 - ASTM D 5182 - DIN 51 354 - CEC L-07-A-71

Datasheet:
15-0249.pdf
Ref.-No:
15-0249 (380/415V)





Congealing Point by Rotating Thermometer (manual)

ASTM D 938 - IP 76 - ISO 2207 (BS 5088 obs.)

Datasheet:
17-0000.pdf

Ref.-No:
17-0000 (230V)
17-0008 (115V)



Melting Point - Cooling Curve (manual)

ASTM D 87 - BS 4695 - FTM 791-1402 - IP 55 - ISO 3841 - JIS K2235

Datasheet:
17-0031.pdf

Ref.-No:
17-0031 (230V)
17-0032 (115V)



Melting Point - Drop (manual)

ASTM D 127 - FTM 791-1401 - IP 133 - JIS K2235

Datasheet:
17-0040.pdf

Ref.-No:
17-0040 (230V)
17-0041 (115V)



Saponification Number (manual)

ASTM D 94-A - DIN 51 559-2 - FTM 791-5401 - IP 136-1 -
ISO 6293 - JIS K2503

Datasheet:
17-0610.pdf

Ref.-No:
17-0610 (230V)
17-0611 (115V)



Neutralization Number by Color-Indicator Titration (manual)

ASTM D 128-18+20 - IP 37 - JIS K2220 (DIN 51 809-1 obs.)

Datasheet:
17-0620.pdf

Ref.-No:
17-0620 (230V)
17-0621 (115V)



Dropping Point (manual)

ASTM D 566 - BS 2877 - DIN 51 801-1 - FTM 791-1421 - IP 132 -
ISO 2176 - JIS K2220

Datasheet:
17-0081.pdf

Ref.-No:
17-0081 (230V)
17-0082 (115V)



Dropping Point by Wide Temperature Range (manual)

ASTM D 2265 - DIN 51 825 - ISO 6299 - JIS K2220

- 6-place aluminum block oven

Datasheet:
17-0151.pdf

Ref.-No:
17-0151 (230V)
17-0152 (115V)



Oil Separation During Storage - ASTM (manual)

ASTM D 1742 - FTM 791-322 - JIS K2220

Datasheet:
17-0111.pdf

Ref.-No:
17-0111 (230V)
17-0112 (115V)

**Oil Separation During Storage - DIN/IP (manual)**

DIN 51 817 - IP 121

Datasheet:
17-0140.pdfRef.-No:
17-0140**Oil Separation by Conical Sieve (manual)**

ASTM D 6184 - FTM 791-321

Datasheet:
17-0131.pdfRef.-No:
17-0131**Oxidation Stability (RSSOT) - PetroOXY (automatic)**

see: Fuel

page 17

Oxidation Stability - OFA (manual)

ASTM D 942 - DIN 51 808 - FTM 791-3453 - VV-L-791e-5314 - IP 142 - JIS K2220

- screw-cap oxidation vessel OFA-2 with test glasses
- digital manometer PA4-OFA
- 4-place bath (liquid "OFA" or dry-heat "OBA/OFA")

Datasheet:
17-0212.pdfRef.-No:
17-0212 (230V)**Flow Pressure by KESTERNICH (manual)**

DIN 51 805

17-0752 for test temperature of -20 °C
17-0754 for test temperature of -35 °CDatasheet:
17-0752.pdfRef.-No:
17-0752 (230V)
17-0754 (230V)**Low-Temperature Torque (automatic)**

ASTM D 1478 (obs.)

Datasheet:
99-0505.pdfRef.-No:
99-0505 (230V)**Evaporation Loss by Wide-Temperature Range (manual)**

ASTM D 2595 - JIS K2220

Datasheet:
17-0390.pdfRef.-No:
17-0390 (230V)**Deleterious Particle Content by Estimation (manual)**

ASTM D 1404

Datasheet:
17-1530.pdfRef.-No:
17-1530



Leakage Tendencies by High Temperature (manual)

ASTM D 3527 - ASTM D 4290 - ASTM D 4950

Datasheet:
99-0106.pdf

Ref.-No:
99-0106 (115/230V)



Corrosion Preventive Properties by Run-In Stand (manual)

ASTM D 1743 - ASTM D 4950

Datasheet:
17-1240.pdf

Ref.-No:
17-1240 (230V)



Corrosion Preventive Properties by EMCOR (manual)

ASTM D 6138 - DIN 51 802 - IP 220 - ISO 11 007

Datasheet:
17-0430.pdf

Ref.-No:
17-0430 (230V)



Resistance to Water Spray (manual)

ASTM D 4049

Datasheet:
17-0465.pdf

Ref.-No:
17-0465 (230V)



Water Washout Characteristics - Static Test (manual)

DIN 51 807-1

- 6-place constant temperature bath

Datasheet:
17-0470.pdf

Ref.-No:
17-0470 (230V)



Mineral Oil & Soap Content by Dialysis (manual)

DIN 51 814

Datasheet:
17-0540.pdf

Ref.-No:
17-0540 (230V)
17-0541



Oil Content (manual)

ASTM D 721 - ASTM D 3235 - IP 158 - FTM 791-5431 - JIS K2235
(DIN 51 571-2, ISO 2908 obs.)

- glass accessories, electric stirrer & vacuum pump
- 3-place bath for dry-ice with cooling coil (17-0550)
- 3-place bath for dry-ice with cooling coil and cryostat (17-0551)

Datasheet:
17-0550.pdf

Ref.-No:
17-0550 (230V)
17-0551 (230V)



Grease Penetration by Cone - PNR 12 (digital)

ASTM D 217 - ASTM D 937 - ASTM D 1321 - ASTM D 1403 -
DAB 10-V.5.8.1 - DIN 51 579 - DIN 51 580 - IP 50 - IP 179 -
IP 310 - IP 376 - ISO 2137 - JIS K2220 - JIS K2235 - NF T60-132 -
SHELL SMS 658 (DIN 51 804-2 obs.)

- cones (Standard, Quarter-Scale, Half-Scale, Micro-cones)
- disk penetrator
- needle penetrator

Datasheet:
18-1120_Grease.pdf
18-1120.pdf

Ref.-No:
18-1120 (100-240V)



Grease Working by Machine (automatic)

ASTM D 217 - FTM 791-313 - IP 50 - ISO 2137 - JIS K2220 -
NF T60-132

- two-place model
- worker with 51-hole or 270-hole worker plate
- worker with screw or bayonet-coupling

Datasheet:
17-0506.pdf

Ref.-No:
17-0506 (230V)



Grease Working by Grease Worker Slave (manual)

ASTM D 217 - IP 50 - ISO 2137 - JIS K2220 - NF T60-132

Datasheet:
17-0500.pdf

Ref.-No:
17-0500



Grease Working by Micro-Grease Worker Quarter-Cone (manual)

ASTM D 1403 - IP 310 - ISO 2137 - JIS K2220 - NF T60-132

Datasheet:
17-0500.pdf

Ref.-No:
18-0134



Grease Working by Micro-Grease Worker SHELL-Cone (manual)

ASTM D 1403 - IP 310 - ISO 2137 - JIS K2220 - NF T60-132

Datasheet:
17-0500.pdf

Ref.-No:
18-0145



Roll-Stability of Lubricating Grease (semi-automatic)

ASTM D 1831 - MIL-G-10924

Datasheet:
17-1600.pdf

Ref.-No:
17-1600 (230V)





Breaking Point by FRAASS - BPA 5 (automatic)

IP 80 - EN 12 593 (DIN 52 012, NF T66-026 obs.) - JIS K2207

- programs for permanent bending, rapid test, etc.
- build-in pelitier cooler
- PC-software BPAcon available

Datasheet:
10-0460.pdf

Ref.-No:
10-0460 (230V)
10-0461 (115V)



Breaking Point by FRAASS (manual)

IP 80 - EN 12 593 (DIN 52 012, NF T66-026 obs.) - JIS K2207

Datasheet:
10-0000.pdf

Ref.-No:
10-0000



Tensile Properties - Force-Ductilometer Method - DDA 3 (automatic)

Elongation: ASTM D 113 - AASHTO T 51 - JIS K2207 - DIN 52013
 Recovery: ASTM D 6084 - EN 13 398 - IP 516
 Force: EN 13 703 - AASHTO T300 - IP 515
 EN 13 589 - IP 520

- bath length of 100 cm or 150 cm available
- 3 programs (semi-automatic for forces < 20 mN; automatic for forces > 20 mN up to 300 N; elastic recovery)
- PC-software DDAcon available

Datasheet:
10-0360.pdf

Ref.-No:
- 100 cm
10-0360 (230V)
10-0361 (115V)

- 150 cm
10-0364 (230V)
10-0365 (115V)



Tensile Properties - Ductilometer Method - DD 3 (semi-automatic)

Elongation: ASTM D 113 - AASHTO T 51 - JIS K2207 - DIN 52013
 Recovery: ASTM D 6084 - EN 13 398 - IP 516

- bath length of 100 cm or 150 cm available
- 3 simultaneous test runs possible
- printer & 3 molds included

Datasheet:
10-0350.pdf

Ref.-No:
- 100 cm
10-0350 (230V)
10-0351 (115V)

- 150 cm
10-0354 (230V)
10-0355 (115V)



Softening Point by Ring & Ball - RKA 5 (automatic)

ASTM D 36 - AASHTO T53 - EN 1427 (DIN 52011, NF T66-008 obs.) -
 IP 58 - JIS K2207
 EN 1871-F
 EN 13 179-1
 ISO 4625-1 - BS 6782-6
 ASTM E 28
 EN 1238

- optional: IP 58 and "Wilhelmi-Method" DIN 1996 pt.15 accessories
- available with electromagnetic ball dispenser and automatic ball application

Datasheet:
10-0800.pdf

Ref.-No:
10-0800 (230V)
10-0801 (115V)



Softening Point by Ring & Ball (manual)

ASTM D 36 - EN 1427 - IP 58 - JIS K2207
 (DIN 52 011, NF T66-008 obs.)
 EN 13 179-1

Datasheet:
10-0281.pdf

Ref.-No:
10-0281 (230V)
10-0284 (115V)



Softening Point by WILHELMI (manual)

DIN 1996-15 - EN 1871-F

Datasheet:
10-0071.pdf

Ref.-No:
10-0071 (230V)



Softening Point by KRAEMER-SARNOW (manual)

DIN 52 025 (DIN 1995 obs.)

Datasheet:
10-0081.pdf

Ref.-No:
10-0081 (230V)
10-0088 (115V)



Effect of Heat and Air - Rolling Thin Film Oven Test Method - RTFOT (semi-automatic)

ASTM D 2872 - Superpave®
EN 12 607-1

Datasheet:
10-0619.pdf

Ref.-No:
10-0619 (230V)



Effect of Heat & Air by Rotating Flask - RFT (manual)

EN 12 607-3 - IP 460-3

Datasheet:
10-0110.pdf

Ref.-No:
10-0110 (230V)



Effect of Heat & Air by Thin-Film Oven - TFOT & Loss on Heating (manual)

ASTM D 1754 - EN 12 607-2 (NF T66-047 obs.) - IP 460-2 - JIS K2207
ASTM D 6 (AASHTO T47 obs.) - EN 13 303 (IP 45, NF T66-011 obs.) - IP 506 - JIS K2207

Datasheet:
10-0117_LOH.pdf
10-0117_TFOT.pdf

Ref.-No:
10-0117 (230V)
10-0118 (115V)



Paraffin Content by Distillation (manual)

EN 12 606-1 (DIN 52 015 obs.) - IP 459-1

Datasheet:
10-0160.pdf

Ref.-No:
10-0160 (230V)



Distillation of Pitch

see: Distillation

page 10

Distillation of Cut-Back Bitumen - ASTM Type (manual)

ASTM D 402 (AASHTO T79 obs.) - EN 13 358 (NF T66-003 obs.) - IP 27

Datasheet:
10-0213.pdf

Ref.-No:
10-0213 (Gas)



Distillation of Cut-Back Bitumen - DIN Type (manual)

DIN 52 024

Datasheet:
10-0173.pdf

Ref.-No:
10-0173 (Gas)



Distillation of Road Tar / Bitumen - ASTM Type (manual)

ASTM D 20 - (ASTM D 246 obs.)

Datasheet:
10-0203.pdf

Ref.-No:
10-0203 (Gas)



Distillation of Road Tar / Bitumen - DIN Type (manual)

DIN 52 027

Datasheet:
10-0152.pdf

Ref.-No:
10-0152 (Gas)



Water Content by Asphalt Distillation - ASTM (manual)

ASTM D 244 - AASHTO T 59 - IP 291 (BS 434 obs.)
- available with flask heating mantle (16-1844+16-1454)

Datasheet:
16-0844.pdf

Ref.-No:
16-0844 (230V)
16-1444 (115V)
16-1844 (230V)
16-1454 (115V)





Water Content by Bitumen Distillation - DIN (manual)

DIN 1996-5
- available with flask heating mantle (16-1845+16-1455)

Datasheet:
16-0845.pdf
Ref.-No:
16-0845 (230V)
16-1445 (115V)
16-1845 (230V)
16-1455 (115V)

page 9

Viscosity by Flow Cups

see: Viscosity



Viscosity, Specific by ENGLER (semi-automatic)

ASTM D 1665 - JIS K2208
IP 212 (DIN 51 560 obs.)

Datasheet:
21-0030.pdf
Ref.-No:
21-0030 (230V)



Viscosity, Apparent by Road & Tar Viscometer (manual)

EN 12 846
- single-unit
- electrically heated
- test cups with different orifices

Datasheet:
10-0312.pdf
Ref.-No:
10-0312 (230V)
10-0313 (115V)



Viscosity by Rotational Viscometer - RV

ASTM D 4402 - AASHTO T 316 (AASHTO TP 48 obs.) -
SHRP 1006 - Superpave®
EN 13 302

Datasheet:
10-0600.pdf
Ref.-No:
10-0600 (230V)



Viscosity by Dynamic Shear Rheometer - DSR

SHRP B-003 - AASHTO T 315 (AASHTO TP 5 obs.)- ASTM P 246 -
SHRP 1007 - Superpave®
EN 12 591 (NF T65-001 obs.)

Datasheet:
10-0770.pdf
Ref.-No:
10-0770 (230V)



Aging by Pressure Aging Vessel - PAV (automatic)

ASTM D 6521 (ASTM PS 36 obs.) - AASHTO R28 (AASHTO PP 1 obs.)
Superpave® - AASHTO M320 (AASHTO MP1 obs.) - SHRP 1001
with CE conformity

Datasheet:
10-0720.pdf
Ref.-No:
10-0720 (230V)
10-0722 (115V)



Stiffness by Bending Beam Rheometer - BBR (automatic)

ASTM D 6648 - AASHTO T 313 (AASHTO TP 1 obs.) - Superpave®
EN 14 771

Datasheet:
10-0700.pdf
Ref.-No:
10-0700 (230V)
10-0701 (115V)
10-0702 (EN, 230 V)

page 12

Flash-Point by TAG Open Cup

see: Flash-Point



Penetration - PNR 12 (digital)

ASTM D 5 - EN 1426 - IP 49 - JIS K2207
- standard needles (precision made, factory or officially certified)
- ball shaped penetrator
- optional cone

Datasheet:
18-1120.pdf
Ref.-No:
18-1120 (100-240V)

page 40

Texture by Rheotex Analyzer - PNR 21

see: Penetration



Thermometers

Thermometers

ASTM / IP

- precision made
- officially certified
- with or without collar

Datasheet:
16-0963.pdf (ASTM)
16-0371.pdf (IP)



Calibration Fluids

Calibration Fluids & Reference Liquids

are available for the following methods:

- Boiling Characteristics
- Coking Tendency by MCRT
- Cloud & Pour Point - CFPP - Freezing Point
- Cold Cranking Simulator
- Density
- Distillation Determination
- Effect of Liquids on Rubber
- Flash-Point Determination
- Fuel Testing
- FZG - Gear Test
- Multi-Test Verification Material (MTVM)
for Bitumen, Gas Oil, Kerosine, Motor Gasoline
- Noack Evaporation Test
- Physical Properties
- Predicting the Borderline Pumping Temperature of Engine Oil
- Salts in Crude Oil
- Shear Stability
- Viscosity
- Viscosity at Low Temperature by BROOKFIELD Viscometer

Datasheet:
26-0702.pdf (Density)
11-0590.pdf (Distillation)
50-9000.pdf (Flash-Point)
15-0256.pdf (Friction Wear)
13-0626.pdf (Fuel)
50-9121.pdf (Miscellaneous)
50-9023.pdf (MTVM, STVM)
16-1080.pdf (Oil)
24-0551.pdf (Viscosity)



Metal Test Specimens

Metal Test Specimens

- Group 40-1000 - Rust Preventives
- Group 40-1500 - Fuels and Fuel Additives
- Group 40-2000 - Solid Film Lubricants
- Group 40-3000 - Antifreeze and Soluble Oils
- Group 40-4000 - Lubricating Oils and Hydraulic Fluids
- Group 40-5000 - Specimens for Special Testing

Datasheet:
40-1001.pdf (Rust)
40-1510.pdf (Fuel)
40-2002.pdf (Solid)
40-3011.pdf (Antifreeze)
40-4010.pdf (Lubricating)
40-5010.pdf (Special)





Texture by Rheotex Analyzer - PNR 21 (digital)

ASTM, IP, ISO, EN, DIN etc. depending on accessories

- control and programming software PNRcon
- Force & Tensile measurement
- additional test-accessories on user's choice

Datasheet:
18-1700.pdf

Ref.-No:
18-1700 (230V)
18-1701 (115V)



Penetration by Penetrometer - PNR 12 (digital)

ASTM, IP, ISO, EN, DIN etc. depending on accessories

- multi-purpose use
- test-accessories on user's choice
- microprocessor controlled with LCD

Datasheet:
18-1120.pdf

Ref.-No:
18-1120 (100-240V)



Penetration by Penetrometer (analog)

ASTM D 5 - ASTM D 217 - ASTM D 1321 - ASTM D 1403 -
ASTM D 5329 - AASHTO T49 - DIN 51 579 - EN 13 880-2 -
IP 49 - IP 50 - IP 179 - IP 310 - IP 376 - ISO 2137 - JIS K2207 -
DAB 10-V.5.8.1 - SHELL SMS 658

Datasheet:
18-0001.pdf

Ref.-No:
18-0001



Penetration Accessories - Test Combinations

sets according standard requirements
incl. needle/cone, plunger, weight, sample container etc.

Datasheet:
18-0132.pdf

Penetration Accessories - Plunger & Weight

18-1160.pdf

Penetration Accessories - Needles

18-0101.pdf

Penetration Accessories - Cones & Discs

18-0260.pdf

Penetration Accessories - Rods & Rams

18-0066.pdf

Penetration Accessories - Sample Jars

18-0065.pdf

Penetration Accessories - Sample Tempering Units

18-2220.pdf

Penetration Accessories - Test Kits

page 35

Penetration Accessories - Grease Worker

Datasheet:
17-0509.pdf
17-0506.pdf



Penetration by Pocket Penetrometer (Field Test)

Datasheet:
99-0015.pdf

Ref.-No:
99-0015

**Oxidation Stability of Food & Cosmetics
- PetroOXY "Stainless Steel" (automatic)**

Correlates with:

ISO 6886 - EN 14112

- with stainless steel test chamber & samle dish

Datasheet:
13-3008.pdf

Ref.-No:

13-3008 (115/230V, EU-plug)

13-3009 (115/230V, EU-plug)



Bending & Breaking Stress - BBS 5 (automatic)

BBS 5 - Standard Temperature Model

Temperature Range: max. approx. -46 to +60°C

Datasheet:
10-0466.pdf

Ref.-No:

10-0466 (230V)

10-0467 (115V)



BBS 5H - High Temperature Model

Temperature Range: max. approx. -46 to +100°C

Datasheet:
10-0466.pdf

Ref.-No:

10-0462 (230V)

10-0463 (115V)



BBS 5C - Climatic Cabinet Model

Temperature Range: depending on climatic cabinet

Datasheet:
10-0466.pdf

Ref.-No:

10-0464 (230V)

10-0465 (115V)



**Water Determination by Titration KARL-FISCHER -
TKF1 (automatic)**

ASTM E 203

Datasheet:
16-3300.pdf

Ref.-No:

16-3300 (90-264V)



**Water Content by Coulometer KARL-FISCHER - CKF1
(automatic)**

ASTM D 1533 - ASTM D 4377 - ASTM D 4928 - ASTM D 6304 -
EN 61 000-3-2 - EN 61 000-3-3 - IP 356 - IP 386 - ISO 10 337 -
ISO 12 937

Datasheet:
16-1880.pdf

Ref.-No:

16-1880 (90-264V)



Water Content by Distillation - DEAN & STARK (manual)

ASTM D 95 - AASHTO T55 - EN 1428 - IP 74 - JIS K2275

- available with flask heating mantle (16-1839+16-1452)

Datasheet:
16-0830.pdf

Ref.-No:

16-0830 (230V)

16-1430 (115V)

16-1839 (230V)

16-1452 (115V)



Water Content by Distillation - DEAN & STARK (extended)

ISO 3733

(also ASTM D 95 - AASHTO T55 - EN 1428 - IP 74 - JIS K2275)

- available with flask heating mantle (16-1850+16-1458)

Datasheet:
16-0850.pdf

Ref.-No:

16-0850 (230V)

16-1450 (115V)

16-1850 (230V)

16-1458 (115V)



Water Content by Distillation - ASTM (manual)

ASTM D 4006 - IP 358 - ISO 9029 - JIS K2275

- available with flask heating mantle (16-1847+16-1457)

Datasheet:
16-0847.pdf

Ref.-No:

16-0847 (230V)

16-1447 (115V)

16-1847 (230V)

16-1457 (115V)





Water Content by Distillation - ASTM (manual)

ASTM D 244 - AASHTO T59 - IP 291
 - available with flask heating mantle (16-1844+16-1454)

Datasheet:
 16-0844.pdf

Ref.-No:
16-0844 (230V)
16-1444 (115V)
16-1844 (230V)
16-1454 (115V)



Water Content by Distillation - DIN (manual)

DIN 1996-5
 - available with flask heating mantle (16-1845+16-1455)

Datasheet:
 16-0845.pdf

Ref.-No:
16-0845 (230V)
16-1445 (115V)
16-1845 (230V)
16-1455 (115V)



Water Content by Distillation - DIN (manual)

DIN 51 432
 - available with flask heating mantle (16-1846+16-1456)

Datasheet:
 16-0846.pdf

Ref.-No:
16-0846 (230V)
16-1446 (115V)
16-1846 (230V)
16-1456 (115V)



Solidification-Point-Apparatus (manual)

ASTM D 852

Datasheet:
 13-1170.pdf

Ref.-No:
13-1170



**Refractive Index & Dispersion
 by Universal ABBE Refractometer (manual)**

DIN 51 423-2

Datasheet:
 13-1520.pdf

Ref.-No:
13-1520 (230V)
13-1521 (115V)



**Refractive Index & Dispersion
 by Modified ABBE Refractometer (manual)**

ASTM D 1218 - ASTM D 1747 - DIN 51 423-2 - ISO 5661 - JIS C2101
 (ASTM D 2159 obs.)

Datasheet:
 13-0285.pdf

Ref.-No:
13-0285 (230V)
13-0300 (115V)



**Refractive Index
 by High Performance Refractometer - ABBEMAT (automatic)**

13-1585 **HP** for substances with a **low refractive index**
 13-1586 **WR** for substances with a **high refractive index**

Datasheet:
 13-1585.pdf

Ref.-No:
13-1585 (100/240V)
13-1586 (100/240V)



**Refractive Index
 by Hand Refractometer - Eclipse (manual) (Field Test)**

ASTM D 5006 - JIS C2101
 for BATTERY (13-1532), FUEL (13-1533) & AVIATION (13-1534) samples

Datasheet:
 13-1532.pdf

Ref.-No:
13-1532
13-1533
13-1534



**Refractive Index
 by Refractometer - RFM 340 (automatic)**

ASTM D 5006 - JIS C2101

Datasheet:
 13-1540.pdf

Ref.-No:
13-1540



Tables of Contents


Title 44 - 47

Order Numbers 48 - 49


Standards 50 - 53



A	page	C	page
Aging by Pressure Aging Vessel - PAV (automatic)	38	Corrosion	
Aging Stability by Baader (manual)	26	- Copper Strip (manual)	14
Air Release Value by Impinger (semi-automatic)	29	- Copper Strip Tarnish Test (manual)	18
Aniline Point		- Copper Strip Test (manual)	18
- Method A (Test Tube) (manual)	15	- Silver Strip Test (manual)	18
- Method B (Thin Film) (manual)	15	- of Engine Coolants (manual)	29
- Method D (U-Tube) (manual)	15	Corrosion Preventing Characteristics	
- Method E (automatic)	15	- Chip & Filterpaper Test (manual)	26
Ash by Muffle Furnace (manual)	26	- HERBERT Corrosion Test (manual)	26
Asphaltenes by Extraction	19, 21	Corrosion Preventive Properties	
Auto-Ignition Temperature (semi-automatic)	13, 15	- EMCOR (manual)	34
Auto-Ignition Temperature - ZPA 3 (automatic)	12, 15	- Run-In Stand (manual)	34
Auto-Ignition Temperature of Solids		D	page
acc. Grewer (semi-automatic)	13, 15	Deleterious Particle Content by Estimation (manual)	33
B	page	Demulsibility Characteristics (manual)	23
Bending & Breaking Stress - BBS 5 (automatic)	41	Demulsibility Characteristics by	
Bending Beam Rheometer - BBR (automatic)	38	HERSCHEL Stirring Method (semi-automatic)	23
Borderline Pumping Temperature (automatic)	26	Demulsification Number (manual)	23
Breaking Point by FRAASS (manual)	36	Density Balance	6, 15
Breaking Point by FRAASS - BPA 5 (automatic)	36	Density or Relative Density	
C	page	- Pressure Hydrometer (manual)	14
Calibration Fluids & Reference Liquids	39	- SCHILLING Effusimeter (manual)	14
Centrifuge	22, 24	- Gases by SCHILLING Effusimeter	6
Cloud and Pour Point		- Light Hydrocarbons by Pressure Hydrometer	6
- CAPP 1 & 4 - Dewar (manual)	25	Depentanization by PONA	15, 21
- CAPP I - Cryostat (semi-automatic)	25	Dielectric Breakdown Voltage (automatic)	29
- CAPP V - Cryostat (semi-automatic)	25	Distillation at Atmospheric Pressure	
Coking Tendency		- ADU 4-DryPoint (automatic)	10
- Carbon Residue	30	- ADU 4 (automatic)	10
- CONRADSON Carbon Residue (manual)	18	- DU 4-Economy (manual)	10
- Micro Carbon Residue - MCR-210 (manual)	18	- DU 4-Professional (manual)	10
- Panel Coking (manual)	30	Distillation of Cut-Back Bitumen	
- RAMSBOTTOM Carbon Residue - RCRT (manual)	18	- ASTM Type (manual)	37
Cold Filter Plugging Point - CFPP / CFPP I / CFPP V (manual)	16	- DIN Type (manual)	37
Color		- Pitch	10, 37
- ASTM-Scale Comparator (manual)	23	- Road Tar / Bitumen - ASTM Type (manual)	37
- Colorimeter - PFX 195/2 (automatic)	23	- Road Tar / Bitumen - DIN Type (manual)	37
- Colorimeter - PFX 880/P (automatic)	23	Dropping Point (manual)	32
- Colorimeter - PFX 995/P (automatic)	24	Dropping Point by Wide Temperature Range (manual)	32
- SAYBOLT-Scale Chromometer	18	Dryness by Freeze Valve (manual)	14
- SAYBOLT-Scale Chromometer (manual)	23		
Compatibility of Elastomer or Rubber with			
Lubricating Fluids & Greases	29		
Conductivity			
- Lab Method - Precision (automatic)	17		
- Portable Method - MLA & DCM (manual) (Field Test)	17		
Congealing Point by Rotating Thermometer (manual)	32		
Contamination by Diaphragm Filtration (manual)	22		



E	page	H	page
Effect of Heat & Air		Heat of Combustion by Calorimeter - C 2000 (automatic)	19
- Rolling Thin Film Oven Test Method - RTFOT (semi-automatic)	37	Hydrogen Sulfide in Liquefied Petroleum (manual)	14
- Rotating Flask - RFT (manual)	37	Hydrolytic Stability by Beverage Bottle (manual)	29
- Thin-Film Oven - TFOT & Loss on Heating (manual)	37	Hydrometer	
Evaporation Loss		- API Specific Gravity - ASTM E 100	5
- NOACK - VP 4000 (automatic)	29	- ASTM Specific Gravity - ASTM E 100	5
- Wide-Temperature Range (manual)	33	- Density - Series L 50 SP	5
- of Lubricants (manual)	29	- Density - Series L 50 SP	5
		- Density - Series M 50 SP	5
F	page	- DIN 12 791 - Series M 100	5
Flash-Point Testing		- Specific Gravity - Series L 50 SP	6
- Sample Changer with PMA 4 (automatic)	11	- Specific Gravity - Series M 50 SP	6
- Closed/Open Cup acc. RAPID Tester "Series 3" (automatic)	12	- (Thermo-), Specific Gravity - ASTM 101H	6
- Closed Cup acc. PENSKEY-MARTENS - PM (manual)	13		
- Closed Cup acc. ABEL (semi-automatic)	12	I	page
- Closed Cup acc. ABEL - ABA 4 (automatic)	11	Insoluble Matters & Total Contamination by Diaphragm Filtration ...	26
- Closed Cup acc. PENSKEY-MARTENS - Army Check (manual)	13		
- Closed Cup acc. PENSKEY-MARTENS - PMA 4 (automatic)	11	L	page
- Closed Cup acc. PENSKEY-MARTENS - PMP 4 (semi-autom.)	12	Lead & Salt by Hot-Extraction - TEL	16, 21
- Closed Cup acc. PENSKEY-MARTENS - PM 4 (semi-automatic)	13	Lead in Gasoline - Iodine Monochloride Method	16, 21
- Closed Cup acc. TAG (semi-automatic)	12, 15	Leakage Tendencies by High Temperature (manual)	34
- Closed Cup acc. TAG - TAG 4 (automatic)	11, 15	Low-Temperature Torque (automatic)	33
- Open Cup acc. TAG (semi-automatic)	12, 15, 38	Lube Oil - VARIO Check (Field Test)	30
Flash- & Fire-Point Testing		Lubrication & Load Capability by BRUGGER (manual)	31
- Open Cup acc. CLEVELAND (manual)	13		
- Open Cup acc. CLEVELAND - CL 5 (semi-automatic)	13	M	page
- Open Cup acc. CLEVELAND Open Cup - CLA 4 (automatic)	11	Melting Point - Cooling Curve (manual)	32
Flock Point (manual)	24	Melting Point - Drop (manual)	32
Flow Pressure by KESTERNICH (manual)	33	Mercaptan by Doctor-Test	19, 21
Fluidity by U-Tube Method (manual)	24	Mercaptan Sulfur by pH-Meter (manual)	21
Fluorescent Indicator Adsorption - FIA	15, 21	Metal Test Specimens	39
Foaming Characteristics - Quadruple (manual)	29	Mineral Oil & Soap Content by Dialysis (manual)	34
Foaming Tendencies of Engine Coolants (manual)	29		
Freezing Point (manual)	16	N	page
Freezing Point (semi-automatic)	16	Neutralization Number	
Freezing Point of Antifreeze Solutions (manual)	16	- Color-Indicator Titration (manual)	32
Friction & Wear by REICHERT M2 (manual)	31	- Color-Indicator Titration (Acid & Base No.) - TAN & TBN	30
		- Potentiometric Titration (Acid & Total Base No.) - TAN & TBN	30
G	page		
Gasoline Diluent of Engine Oil by Distillation (manual)	26	O	page
Grease Penetration by Cone - PNR 10 (digital)	35	Oil Content (manual)	34
Grease Penetration by Cone - PNR 12 (digital)	35	Oil Content of Water Mixed Cooling Lubricants by Acidification & Saponification (manual)	25
Grease Working		Oil Separation by Conical Sieve (manual)	33
- Grease Worker Slave (manual)	35	Oil Separation During Storage - ASTM (manual)	32
- Machine (automatic)	35	Oil Separation During Storage - DIN/IP (manual)	33
- Micro-Grease Worker Quarter-Cone (manual)	35		
- Micro-Grease Worker SHELL-Cone (manual)	35		
GUM by Jet Evaporation (manual)	15		



O	page	R	page
Oxidation Characteristics		Refractive Index & Dispersion by	
- TOST/TOO Dry-Heat - ASTM (manual)	27	- Modified ABBE Refractometer	16, 42
- TOST/TOO Dry-Heat - DIN (manual)	27	- Universal ABBE Refractometer	16, 42
- TOST/TOO Liquid- ASTM D 2893 (manual)	27	Refractive Index	
- TOST/TOO Liquid - ASTM D 943 (manual)	27	Hand Refractometer - Eclipse	16, 42
Oxidation Stability & Corrosiveness (manual)	27	High Performance Refractometer - ABBEMAT	16, 42
Oxidation Stability		Refractometer	30, 42
- BAM 200, BAM 250	27, 28	Refractometer - RFM 340 (automatic)	16, 42
- BAM Dry-Heat (manual)	28	Resistance to Water Spray (manual)	34
- BAM Liquid (manual)	27	Resistance by Philipp Test (manual)	24
- Dry-Heat Bath Safety Test Arrangement - OBA 2/T	18	Roll-Stability of Lubricating Grease (semi-automatic)	35
- Dry-Heat Bath Test Arrangement - OBA 1/T	17	Rust Preventing by Chip & Filter (manual)	28
- Liquid Bath Safety Test Arrangement - OBA 2	18	Rust Preventing Characteristics - TOR (semi-automatic)	28
- Liquid Bath Test Arrangement - OBA 1	17	Rust Protection by Humidity Cabinet (manual)	28
- OBA (semi-automatic)	17		
- OFA (manual)	33	S	page
- RSSOT - PetroOXY (automatic)	17, 27, 33	Salts by Conductivity (electrometric)	
- of Food/Cosmetics - PetroOXY (automatic)	41	- Salinometer Portable (manual)	19, 21
- Rotating Bomb - RBOT (automatic)	28	Salts by Conductivity (electrometric)	
- Rotating Bomb - RBOT (semi-automatic)	28	- Salinometer SCO 1 (manual)	19, 21
		Sampling	
P	page	- Cylinders - double valve	4
Paraffin Content by Distillation (manual)	37	- Devices	4
Particulates & Filter Membrane Color Ratings in		- Tubes - closed	4
Aviation Fuel In-Line Sampler - AFCON	19, 21	- Tubes - open	4
Particulates & Free Water	19, 20	- Tubes - open (Grease)	4
Particulate Contamination in Aviation Fuel by Laboratory Filtration	19	Saponification Number (manual)	32
PCB-Content - PCB Check (Field Test)	30	Sediment, Total (manual)	22
Penetration - PNR 12(digital)	38	Sediment by Extraction	19, 22
Penetration Accessories		Shear Stability by Diesel Injector Rig (manual)	25
- Cones & Discs	40	Shear Stability of Polymer-Containing Oils and Hydraulic Fluids	
- Grease Worker	40	- Sonic (semi-automatic)	25
- Needles	40	Smoke Point by Lamp Method (manual)	18
- Plunger & Weight	40	Softening Point	
- Rods & Rams	40	- KRAEMER-SARNOW (manual)	36
- Sample Jars	40	- Ring & Ball (manual)	36
- Sample Tempering Units	40	- Ring & Ball - RKA 5 (automatic)	36
- Test Combinations	40	- WILHELMI (manual)	36
- Test Kits	40	Solidification-Point-Apparatus (manual)	42
Penetration by Penetrometer		Solidity & Turbidity Point (manual)	24
- analog	40	Stability	
- PNR 12 (digital)	40	- Oil / Layer Separation (manual)	25
- PNR 21 (Rheotex)	40	- Thermal for Jet Fuel (manual)	19
- Pocket Penetrometer (Field Test)	40	Stiffness by Bending Beam Rheometer - BBR (automatic)	38
pH-Value by pH-Meter (manual)	26	Sulfonation Number by Shaking Machine (manual)	20
Pycnometer			
- GAY-LUSSAC	6		
- HUBBARD	6		
- LIPKIN	6		

**S** page

Sulfur & Chlorine	
- Active	20, 26
- GROTE-KREKELER (manual)	20
- High Temperature (manual)	20
- WICKBOLD Combustion	19, 20
- Raney-Nickel (manual)	20
- SCHÖNIGER (Thorin-Titration) (manual)	20

T page

Tank Sampling Bottles	4
Tank Sampling by Scoop-Thermometer	4
Tempering Baths	
- Density, API & Specific Gravity	5
- Viscosity - Kinematic Capillary Methods	7
- Viscositytubes - Ambient Temperature	7
- Viscositytubes - High Temperature	7
Tensile Properties - Force-Ductilometer Method - DDA 3 (automatic)	36
Tensiometer	
- Processor (automatic)	24
- Scale (manual)	24
Texture by Rheotex Analyzer - PNR 21	38, 40
Thermometers	39
Total Dirt by Refrigerated Centrifuge (semi-automatic)	24
Total Sediment	19
Tensile Properties - Ductilometer Method - DD 3 (semi-autom.)	36

U page

Unulfonated Residue of Plant-Spray-Oil by Combined Shaking Machine (manual)	20
---	----

V page


Vapor Pressure (manual)	14
Vapor Pressure Analyzer (automatic)	16
Vapor Pressure Analyzer - PetroVap II (automatic)	16
Vapor Pressure by REID (digital)	15
Viscometer Tubes	
- BS/IP - BS/IP/SL(S) (Suspended Level)	8
- BS/IP U-Tube - BS/IP/RF (Reverse Flow)	8
- CANNON-FENSKE Opaque - CFOV (Reverse Flow)	8
- CANNON-FENSKE Routine - CFRV (Ostwald)	8
- CANNON-UBBELOHDE Semi-Micro - CUSMV (Suspended Level)	7
- Micro-OSTWALD - MOV (Ostwald)	8
- UBBELOHDE ASTM-Type - UAV (Suspended Level)	7
- UBBELOHDE DIN-Type - UDV (Suspended Level)	7
- UBBELOHDE Micro Viscometers - UMV (Suspended Level)	7
- Vacuum - CMVV, AIVV, MKVV	8
- ZEITFUCHS Cross-Arm - ZCAV (Reverse Flow)	8
- UBBELOHDE ASTM/DIN Capillary - ViscoClock (semi-autom.)	7
- UBBELOHDE Capillary - AKV (automatic)	7

V page

Viscometer	
- Apparent by Road & Tar Viscometer (manual)	38
- Brookfield (semi-automatic & automatic)	25
- Comparator - Viscosity Check	9
- Comparator - Viscosity Check (Field Test)	30
- Constant Shear Stress Viscometer	26
- Dynamic Shear Rheometer - DSR	38
- Falling Ball - HÖPPLER (manual)	9
- Falling Ball - JUNG Check	9
- Falling Ball - JUNG Check (Field Test)	30
- Flow Cups - FORD, ISO (manual)	9, 38
- High Temperature & High Shear Rate - HTHS	25
- Low Temperature	25, 26
- Mini-Rotary Viscometer - MRV	26
- Oil Standards	8
- REDWOOD (manual)	9
- Rotational Viscometer - RV	38
- SAYBOLT (manual)	9
- Specific by ENGLER	8
- Specific by ENGLER (semi-automatic)	38
- Tapered-Plug Viscometer - HTHS (semi-automatic)	25
Volatility & Residue by Propane Residue (Mercury Freeze Method) (manual)	14

W page

Water & Sediment by Centrifuge (semi-automatic)	19, 22
Water Content	
- Asphalt Distillation - ASTM (manual)	30, 37
- Bitumen Distillation - DIN (manual)	30, 38
- Distillation - ASTM (manual)	30, 41, 42
- Distillation - DEAN & STARK (extended)	30, 41
- Distillation - DEAN & STARK (manual)	30, 41
- Distillation - DIN (manual)	30, 42
- Coulorimeter - KARL-FISCHER - CKF1 (automatic)	30, 41
- Titration - KARL-FISCHER - TKF1 (automatic)	41
Water in Oil - WIO Check (Field Test)	30
Water Reaction by Mixing Cylinder (manual)	17
Water Separation Characteristics by Micro-Separometer - MSEP (Field Test)	17
Water Vapor Content by Dew-Point (manual)	14
Water Washout Characteristics - Static Test (manual)	34
Wear Characteristics	
- FALEX Pin & Vee Block (manual)	31
- SHELL - Four-Ball Machine (manual)	31
- TIMKEN-Test (manual)	31
- FZG Gear Oil Rig (manual)	31



10-	page	12-	page	13-	page	14-	page
10-0000	36	12-0501 (115/230V)	11	13-0870 (230V)	18	14-0004	14
10-0071 (230V)	36	12-0503 (115/230V)	11	13-0871 (115V)	18	14-0015	14
10-0081 (230V)	36	12-0540 (115/230V)	11	13-0920 (230V)	20	14-0047 (115V)	14
10-0088 (115V)	36	12-0542 (115/230V)	11	13-0929 (115V)	15	14-0050 (230V)	14
10-0110 (230V)	37	12-0595 (115/230V)	12	13-0930 (230V)	15	14-0060	14
10-0117 (230V)	37	12-0601 (Gas)	13	13-1000 (230V)	18	14-0090	14
10-0118 (115V)	37	12-0701 (Gas)	13	13-1019 (115V)	18	14-0100	14
10-0152 (Gas)	37	12-1000 (230V)	13	13-1020 (230V)	18	14-0110 (230V)	14
10-0160 (230V)	37	12-1001 (115V)	13	13-1028 (230V)	18	14-0116 (°F)	6
10-0173 (Gas)	37	12-1240 (115/230V)	12	13-1029 (115V)	18	14-0121 (°C)	6
10-0179 (230V)	10	12-1241 (115/230V)	12	13-1030 (230V)	15	14-0230 (230V)	14
10-0203 (Gas)	37	12-1260 (115/230V)	12	13-1031 (115V)	15	14-0231 (115V)	14
10-0213 (Gas)	37	12-1270 (230 V)	12	13-1160 (230V)	20		
10-0281 (230V)	36	12-1520 (230V)	12	13-1161 (115V)	20		
10-0284 (115V)	36	12-1521 (110V)	12	13-1170	42		
10-0312 (230V)	38	12-1540 (230V)	12	13-1410 (Battery)	17	15-	page
10-0313 (115V)	38	12-1650 (230V)	13	13-1470 (230V)	16	15-0035 (230V)	31
10-0350 (230V)	36	12-1767 (115/230V)	11	13-1490 (230V)	21	15-0070 (380V)	31
10-0351 (115V) - 100 cm	36	12-1770 (115/230V)	11	13-1491 (115V)	21	15-0249 (380/415V)	31
10-0354 (230V)	36	12-1802 (230V)	13	13-1520 (230V)	42	15-0450 (230V)	31
10-0355 (115V) - 150 cm	36	12-1823 (115V)	12	13-1521 (115V)	42	15-0602 (230V)	31
10-0360 (230V)	36	12-1840 (230V)	13	13-1532	42	15-0603 (115V)	31
10-0361 (115V) - 100 cm	36	12-2730 (230V)	13	13-1533	42	15-1000 (380V)	31
10-0364 (230V)	36	12-2731 (115V)	13	13-1534	42		
10-0365 (115V) - 150 cm	36	12-2770 (with EU-plug)	11	13-1540	42		
10-0460 (230V)	36	12-2771 (with US-plug)	11	13-1560 (230V)	21	16-	page
10-0461 (115V)	36	12-2772 (with EU-plug)	11	13-1585 (100/240V)	42	16-0001 (230V)	26
10-0462 (230V)	41	12-2773 (with US-plug)	11	13-1586 (100/240V)	42	16-0070 (230V)	25
10-0463 (115V)	41	12-2776 (with EU-plug)	11	13-1650	30	16-0072 (230V)	25
10-0464 (230V)	41	12-2777 (with US-plug)	11	13-1660 (230V)	17	16-0074 (230V)	25
10-0465 (115V)	41	12-2777 (with EU-plug)	11	13-1661 (115V)	17	16-0082 (230V)	25
10-0466 (230V)	41	12-2778 (with US-plug)	11	13-1662 (230V)	17	16-0084 (230V)	20
10-0467 (115V)	41	12-5700 (110/230V)	11	13-1663 (115V)	17	16-0085 (115V)	20
10-0600 (230V)	38			13-1664 (230V)	18	16-0140 (230V)	23
10-0619 (230V)	37			13-1665 (115V)	18	16-0141 (115V)	23
10-0700 (230V)	38			13-1666 (230V)	18	16-0179 (115V)	23
10-0701 (115V)	38			13-1667 (115V)	18	16-0180 (230V)	23
10-0702 (230V)	38			13-1720	20	16-0210 (115/230V)	23
10-0720 (115/230V)	38			13-1880 (230V)	19	16-0249	28
10-0770 (230V)	38			13-1980	21	16-0252 (230V)	24
10-0800 (230V)	36			13-2140 (230V)	22	16-0253 (115V)	24
10-0801 (115V)	36			13-2141 (120V)	22	16-0281 (230V)	26
				13-2160 (230V)	22	16-0282 (115V)	26
				13-2161 (120V)	22	16-0286 (230V)	26
				13-2300 (230V)	22	16-0287 (115V)	26
				13-2301 (115V)	22	16-0295	26
				13-2340	21	16-0301	24
				13-2350	21	16-0320 (380V)	25
				13-2414 (115V)	16	16-0352	25
				13-2418 (230V)	16	16-0389 (115V)	24
				13-2419 (230V)	16	16-0390 (230V)	24
				13-2460 (230V)	16	16-0401 (230V)	24
				13-2600 (105-250V)	16	16-0402 (115V)	24
				13-2620 (115/230V)	16	16-0431	21
				13-2900 (230V)	18	16-0458 (115V)	27
				13-3002 (EU-plug) ...	17, 27, 33	16-0581 (230V)	27
				13-3003 (US-plug) ...	17, 27, 33	16-0583 (230V)	27
				13-3006 (EU-plug) ...	17, 27, 33	16-0584 (115V)	27
				13-3007 (US-plug) ...	17, 27, 33	16-0683 (230V)	26
				13-3008 (EU-plug)	41	16-0689 (115V)	29
				13-3009 (US-plug)	41	16-0690 (230V)	29
				13-3180 (Battery)	17	16-0724 (230V)	28
						16-0725 (115V)	28
						16-0749 (115V)	26
						16-0750 (230V)	26
						16-0771 (230V)	29
						16-0772 (115V)	29
						16-0830 (230V)	41
						16-0844 (230V)	37, 42
						16-0845 (230V)	38, 42
						16-0846 (230V)	42
						16-0847 (230V)	41

11-	page	13-	page
11-0540 (230V)	10	13-0001 (230V)	15
11-0541 (115V)	10	13-0035 (230V)	15
11-0542 (230V)	10	13-0036 (115V)	15
11-0543 (115V)	10	13-0100 (DIN)	15
11-0900 (230V)	10	13-0101 (ASTM)	15
11-0901 (115V)	10	13-0102 (230V)	15
11-0902 (230V)	10	13-0103 (115V)	15
11-0903 (115V)	10	13-0106 (230V)	15
11-0904 (230V)	10	13-0107 (115V)	15
11-0905 (115V)	10	13-0169 (115V)	21
11-0906 (230V)	10	13-0170 (230V)	21
11-0907 (115V)	10	13-0195 (Battery)	17
11-0910 (230V)	10	13-0201 (<180 bar)	15
11-0911 (115V)	10	13-0202 (>180 bar)	15
11-0912 (230V)	10	13-0272	19
11-0913 (115V)	10	13-0285 (230V)	42
11-1540 (230V)	10	13-0300 (115V)	42
11-1541 (115V)	10	13-0305 (220V)	21
11-1542 (230V)	10	13-0321	22
11-1543 (115V)	10	13-0391	17
11-5580 (230V)	10	13-0401	17
11-5581 (115V)	10	13-0420	16
11-5582 (230V)	10	13-0461	16
11-5583 (115V)	10	13-0475 (230V)	20
		13-0497 (230V)	21
		13-0670 (230V)	20
		13-0671 (115V)	20
		13-0691	18
		13-0699 (115V)	18
		13-0700 (230V)	18
		13-0756 (230V)	20
		13-0803 (230V)	20
		13-0841 (Gas)	18
		13-0842 (Gas)	18
		13-0843 (Gas)	18



16-	page	16-	page	18-	page	24-	page
16-0850 (230V)	41	16-2360 (230V)	25	18-0001	40	24-0025	9
16-0955 (230V)	29	16-2363 (115V)	25	18-0134	35	24-0501 (Secondary)	8
16-0970	22	16-2372 (230V)	25	18-0145	35	24-0551 (Primary)	8
16-1152 (230V)	25	16-2373 (115V)	25	18-1120 (100-240V)	35, 38, 40	24-0601	8
16-1200 (230V)	29	16-2410 (230V)	25	18-1700 (230V)	40		
16-1201 (115V)	29	16-2430 (220V)	26	18-1701 (115V)	40		
16-1230 (230V)	23	16-2431 (220V)	26			25-	page
16-1317 (230V)	26	16-2432 (110V)	26	19-	page	25-0716 (230V)	27
16-1383 (3x 220V)	25	16-2500 (220V)	25	19-0001	4		
16-1416 (230V)	26	16-2501 (220V)	25	19-0002	4		
16-1417 (115V)	26	16-2502 (110V)	25	19-0003	4		
16-1418 (230V)	26	16-2503 (220V)	25	19-0005	4	26-	page
16-1419 (115V)	26	16-2504 (220V)	25	19-0011	4	26-0010 (230V)	5
16-1420 (230V)	27	16-2740 (60kV)	29	19-0012	4	26-0011 (115V)	5
16-1421 (115V)	27	16-2745 (80kV)	29	19-0032 (BACON)	4	26-0017 (230V)	5
16-1430 (115V)	41	16-2750 (100kV)	29	19-0040 (Can)	4	26-0018 (115V)	5
16-1444 (115V)	37, 42	16-3250 (230V)	28	19-0041 (Glas)	4	26-0019 (115/230V)	6
16-1445 (115V)	38, 42	16-3300 (90-264V)	41	19-0051 (Auger)	4	26-0400	5
16-1446 (115V)	42			19-0056 (ASTM)	4	26-0405	5
16-1447 (115V)	41	17-	page	19-0058 (DIN)	4	26-0500	5
16-1450 (115V)	41	17-0000 (230V)	32	19-0070	4	26-0509	5
16-1452 (115V)	41	17-0008 (115V)	32	19-0071	4	26-0600	5
16-1454 (115V)	37, 42	17-0031 (230V)	32	19-0150	4	26-0608	5
16-1455 (115V)	38, 42	17-0032 (115V)	32	19-0151	4	26-0712	6
16-1456 (115V)	42	17-0040 (230V)	32	19-0152	4	26-0721	6
16-1457 (115V)	41	17-0041 (115V)	32	19-0153	4	26-0740 (IP)	6
16-1458 (115V)	41	17-0081 (230V)	32	19-0154	4	26-0741 (ISO)	6
16-1470 (230V)	29	17-0082 (115V)	32	19-0155	4		
16-1471 (115V)	29	17-0111 (230V)	32	19-0156	4		
16-1472 (230V)	29	17-0112 (115V)	32			40-	page
16-1473 (115V)	29	17-0131	33	20-	page	40-0120	30
16-1500 (230V)	27	17-0140	33	20-0000 (230V)	7	40-0140	30
16-1501 (115V)	27	17-0151 (230V)	32	20-0001 (115V)	7	40-0150	30
16-1502 (230V)	27	17-0152 (115V)	32	20-0150 (230V)	7	40-0160	30
16-1503 (115V)	27	17-0212 (230V)	33	20-0151 (115V)	7	40-1001	39
16-1519 (230V)	28	17-0390 (230V)	33			40-1510	39
16-1520 (115V)	28	17-0430 (230V)	34	21-	page	40-2002	39
16-1522 (230V)	28	17-0465 (230V)	34	21-0030 (230V)	8, 38	40-3011	39
16-1523 (115V)	28	17-0470 (230V)	34	21-0300 (230V)	9	40-4010	39
16-1566 (230V)	28	17-0500	35	21-0335 (230V)	9	40-5010	39
16-1568 (115V)	28	17-0506 (230V)	35	21-0501 (230V)	9		
16-1569 (230V)	28	17-0540 (230V)	34	21-0502 (115V)	9	99-	page
16-1570 (115V)	28	17-0541	34	21-0502 (115V)	9	99-0015	40
16-1660 (230V)	25	17-0550 (230V)	34	21-0750 (DIN)	9	99-0106 (115/230V)	34
16-1750 (115/230V)	24	17-0551 (230V)	34	21-0761 (ASTM)	9	99-0117 (230V)	29
16-1755	24	17-0600	30	21-0766 (ISO)	9	99-0300 (220/240V)	30
16-1839 (230V)	41	17-0610 (230V)	32			99-0505 (230V)	33
16-1844 (230V)	37, 42	17-0611 (115V)	32	22-	page	99-1031	5
16-1845 (230V)	38, 42	17-0620 (230V)	32	22-0000 (manual)	8	99-1040	5
16-1846 (230V)	42	17-0621 (115V)	32	22-0050 (manual)	8	99-1131	5
16-1847 (230V)	41	17-0695 (230V)	30	22-0100 (manual)	7	99-1140	5
16-1850 (230V)	41	17-0696 (115V)	30	22-0120 (manual)	7	99-1141	5
16-1870 (230V)	24	17-0752 (230V)	33	22-0175 (230V)	7	99-1150	5
16-1880 (90-264V)	41	17-0754 (230V)	33	22-0176 (115V)	7	99-1151	6
16-1960	19	17-0780 (230V)	29	22-0200 (manual)	8	99-1160	6
16-1980 (230V)	28	17-1240 (230V)	34	22-0230 (manual)	8	99-1161	6
16-1990 (230V)	28	17-1530	33	22-0300 (automatic)	7	99-1170	6
16-1991 (115V)	28	17-1600 (230V)	35	22-0320 (automatic)	7		
16-2010 (230V)	29	17-1800 (110-240V)	30	22-0360 (manual)	7		
16-2011 (115V)	29			22-0390 (manual)	7		
16-2100 (230V)	23			22-0395 (automatic)	7		
16-2101 (115V)	23			22-0400 (manual)	8		
16-2175 (115/230V)	23			22-0405 (automatic)	8		
16-2176 (115/230V)	23			22-0500	8		
16-2190 (115/230V)	24			22-0520	8		
16-2191 (115/230V)	24			22-0530 (manual)	8		
16-2194 (115/230V)	23			22-1000 (230V)	7		
16-2348 (230V)	25			22-1001 (115V)	7		
16-2349 (115V)	25						




American Standards

AASHTO	page	ASTM	page	ASTM	page	ASTM	page
AASHTO MP 1 (obs.)	38	ASTM D 524	18	ASTM D 1935 (obs.)	23	ASTM D 4206	12
AASHTO M 320	38	ASTM D 525	17	ASTM D 1989 (obs.)	19	ASTM D 4289	29
		ASTM D 526 (obs.)	21	ASTM D 2001	15	ASTM D 4290	34
AASHTO PP 1 (obs.)	38	ASTM D 566	32	ASTM D 2111	5	ASTM D 4308	17
		ASTM D 611-A+B+D+E	15	ASTM D 2112	28	ASTM D 4310	27
AASHTO R28	38	ASTM D 664	30	ASTM D 2158	14	ASTM D 4377	41
		ASTM D 665	28	ASTM D 2159 (obs.)	42	ASTM D 4422	26
AASHTO T47 (obs.)	37	ASTM D 721	34	ASTM D 2162	8	ASTM D 4530	18
AASHTO T48	11, 13	ASTM D 849	18	ASTM D 2170	7, 8	ASTM D 4636	27
AASHTO T49	40	ASTM D 850	10	ASTM D 2171	8	ASTM D 4684	26
AASHTO T51	36	ASTM D 852	42	ASTM D 2265	32	ASTM D 4739 (alternative)	30
AASHTO T53	36	ASTM D 873	17	ASTM D 2266	31	ASTM D 4741	25
AASHTO T55	22, 41	ASTM D 874	26	ASTM D 2272	17, 28	ASTM D 4742	28
AASHTO T59	37, 42	ASTM D 877	29	ASTM D 2273	22	ASTM D 4809	19
AASHTO T72	9	ASTM D 892 (seq. I-III)	29	ASTM D 2274	17, 27	ASTM D 4860	17
AASHTO T73	11 - 13	ASTM D 893	22	ASTM D 2276	21	ASTM D 4870	22
AASHTO T79 (obs.)	37	ASTM D 923	4	ASTM D 2384	20	ASTM D 4928	41
AASHTO T172	11 - 13	ASTM D 937	35	ASTM D 2386	16	ASTM D 4950	34
AASHTO T300	36	ASTM D 938	32	ASTM D 2420	14	ASTM D 4952	21
AASHTO T313	38	ASTM D 942	17, 33	ASTM D 2440	27	ASTM D 4953	15, 16
AASHTO T315	38	ASTM D 943	17, 27	ASTM D 2500	25	ASTM D 4998	31
AASHTO T316	38	ASTM D 971	24	ASTM D 2509	31	ASTM D 5006	42
		ASTM D 972	29	ASTM D 2533	16	ASTM D 5125	9
AASHTO TP1 (obs.)	38	ASTM D 974	30	ASTM D 2547 (obs.)	21	ASTM D 5182	31
AASHTO TP5 (obs.)	38	ASTM D 1026	26	ASTM D 2569	10	ASTM D 5183	31
AASHTO TP48 (obs.)	38	ASTM D 1078	10	ASTM D 2595	33	ASTM D 5190	16
		ASTM D 1094	17	ASTM D 2596	31	ASTM D 5191	16
		ASTM D 1121	26	ASTM D 2603	25	ASTM D 5304	17
		ASTM D 1142	14	ASTM D 2619	29	ASTM D 5329	40
		ASTM D 1177	16	ASTM D 2624	17	ASTM D 5452	19
ASTM	page	ASTM D 1200	9	ASTM D 2625	31	ASTM D 5468	19
ASTM D 5	38, 40	ASTM D 1209	23	ASTM D 2670	31	ASTM D 5482	16
ASTM D 6	37	ASTM D 1218	42	ASTM D 2709	22	ASTM D 5534	28
ASTM D 20	37	ASTM D 1265	4	ASTM D 2711	23	ASTM D 5621	25
ASTM D 36	36	ASTM D 1267	14	ASTM D 2713	14	ASTM D 5800-B	29
ASTM D 56	11, 12	ASTM D 1287	26	ASTM D 2782	31	ASTM D 5853	25
ASTM D 70	5, 6	ASTM D 1298	5	ASTM D 2783	31	ASTM D 6045	23, 24
ASTM D 71	5, 6	ASTM D 1310	12	ASTM D 2784	20	ASTM D 6082 (seq. IV)	29
ASTM D 86	10	ASTM D 1319	15	ASTM D 2872	37	ASTM D 6084	36
ASTM D 87	32	ASTM D 1321	35, 40	ASTM D 2878	29	ASTM D 6138	34
ASTM D 88	9	ASTM D 1322	18	ASTM D 2896	30	ASTM D 6166	23, 24
ASTM D 91	22	ASTM D 1331	24	ASTM D 2983	25	ASTM D 6184	33
ASTM D 92	11, 13	ASTM D 1384	29	ASTM D 3143	12	ASTM D 6278	25
ASTM D 93-A+B	11, 12, 13	ASTM D 1401	23	ASTM D 3227	20	ASTM D 6304	41
ASTM D 94-A	32	ASTM D 1403	35, 40	ASTM D 3230	21	ASTM D 6371	16
ASTM D 95	41	ASTM D 1404	33	ASTM D 3233	31	ASTM D 6378	16
ASTM D 96 (obs.)	22	ASTM D 1417	24	ASTM D 3235	34	ASTM D 6521	38
ASTM D 97	25	ASTM D 1478 (obs.)	33	ASTM D 3241	19	ASTM D 6560	21
ASTM D 113	36	ASTM D 1481	5	ASTM D 3278-A+B	12	ASTM D 6648	38
ASTM D 127	32	ASTM D 1500	23, 24	ASTM D 3300	26	ASTM D 6751	11, 12, 13
ASTM D 128-18+20	32	ASTM D 1533	41	ASTM D 3341	21	ASTM D 6821	26
ASTM D 130	18	ASTM D 1544	23	ASTM D 3427	29	ASTM D 6896	26
ASTM D 156	23, 24	ASTM D 1552	20	ASTM D 3527	34	ASTM D 7109	25
ASTM D 189	18	ASTM D 1590 (obs.)	24	ASTM D 3603	28	ASTM D 7525	17
ASTM D 217	35, 40	ASTM D 1657	6, 14	ASTM D 3828-A+B	12	ASTM D 7545	17
ASTM D 240	19	ASTM D 1661 (obs.)	19	ASTM D 3829	26	ASTM D 7724	17
ASTM D 244	37, 42	ASTM D 1662	20	ASTM D 3934	11		
ASTM D 246 (obs.)	37	ASTM D 1665	38	ASTM D 3941	11	ASTM E 28	36
ASTM D 270 (obs.)	4	ASTM D 1742	32	ASTM D 3948	17	ASTM E 100	5
ASTM D 287	5	ASTM D 1743	34	ASTM D 4006	41	ASTM E 102	9
ASTM D 322	26	ASTM D 1747	42	ASTM D 4007	22	ASTM E 203	41
ASTM D 323	16	ASTM D 1748	28	ASTM D 4048	18	ASTM E 308 (obs.)	23, 24
ASTM D 323-A-B-D	15	ASTM D 1754	37	ASTM D 4049	34	ASTM E 502 (obs.)	12
ASTM D 381	15	ASTM D 1796	22	ASTM D 4052	8	ASTM E 659	13
ASTM D 402	37	ASTM D 1816	29	ASTM D 4057	4		
ASTM D 445	7, 8	ASTM D 1831	35	ASTM D 4172	31	ASTM P 246	38
ASTM D 446	7, 8	ASTM D 1837	14	ASTM D 4176	21		
ASTM D 473	22	ASTM D 1838	14	ASTM D 4177	4	ASTM PS 36 (obs.)	38
ASTM D 482	26	ASTM D 1881	29				
ASTM D 483	20						



British Standards

FTM	page	BS	page	IP	page	IP	page
FTM 141-4293	11, 12, 13	BS 148	26	IP 2-A+B+D+E	15	IP 196	23, 24
FTM 791-101	23	BS 434 (obs.)	37	IP 4	26	IP 212	38
FTM 791-102 (obs.)	23, 24	BS 684-1.17	11, 12, 13	IP 13	18	IP 216	21
FTM 791-201	25	BS 2634 (obs.)	30	IP 14	18	IP 219	25
FTM 791-304	9	BS 2877	32	IP 15	25	IP 220	34
FTM 791-305	7, 8	BS 3900-A11	12	IP 16	16	IP 222	8
FTM 791-313	35	BS 4347 (obs.)	17	IP 17-A+B	24	IP 227 (obs.)	18
FTM 791-321	33	BS 4348 (obs.)	15	IP 19	23	IP 229	17, 28
FTM 791-322	32	BS 4351 (obs.)	18	IP 23	26	IP 235	14
FTM 791-351	29	BS 4380 (obs.)	18	IP 27	37	IP 239	31
FTM 791-380 (obs.)	23	BS 4382 (obs.)	22	IP 30	21	IP 240	31
FTM 791-1040	20	BS 4385 (obs.)	22	IP 34-A+B	11, 12, 13	IP 242	20
FTM 791-1101	11, 12	BS 4388 (obs.)	27	IP 35 (optional)	13	IP 243	20
FTM 791-1102	11, 12, 13	BS 4451 (obs.)	18	IP 36	11, 13	IP 248 (obs.)	21
FTM 791-1103	11, 13	BS 4456 (obs.)	17	IP 37	32	IP 263 (similar)	25
FTM 791-1401	32	BS 4457 (obs.)	30	IP 40	17	IP 265	21
FTM 791-1402	32	BS 4458 (obs.)	25	IP 45 (obs.)	37	IP 267-B	25
FTM 791-1421	32	BS 4695	32	IP 48	27, 28	IP 276	30
FTM 791-2107	18	BS 4699 (obs.)	5, 6	IP 49	38, 40	IP 280	27
FTM 791-2506	19	BS 4714 (obs.)	5	IP 50	35, 40	IP 287	28
FTM 791-3002	22	BS 5088 (obs.)	32	IP 51 (obs.)	4	IP 291	42
FTM 791-3008	21	BS 6782-6	36	IP 53	22	IP 294	25
FTM 791-3302	15	BS 7120	28	IP 55	32	IP 295	29
FTM 791-3352	17			IP 57	18	IP 306	27
FTM 791-3354	17			IP 58	36	IP 307	27
FTM 791-3453	33			IP 59-A+B (obs.)	6	IP 309	16
FTM 791-3462	30			IP 59-C (obs.)	14	IP 310	35, 40
FTM 791-3601	15			IP 69	16	IP 313	29
FTM 791-3703	15			IP 69-A-B-D	15	IP 317	14
FTM 791-4011	28			IP 70 (obs.)	9	IP 319	7, 8
FTM 791-5001	18			IP 71-1+2	7, 8	IP 323	19
FTM 791-5002	18			IP 74	41	IP 326	31
FTM 791-5105	30			IP 75 (obs.)	22	IP 335	27
FTM 791-5106	30			IP 76	32	IP 342	20
FTM 791-5308	27			IP 77	21	IP 356	41
FTM 791-5315	28			IP 80	36	IP 358	41
FTM 791-5325	18			IP 96 (obs.)	21	IP 359 (obs.)	22
FTM 791-5401	32			IP 112	18	IP 375	22
FTM 791-5421	26			IP 121	33	IP 376	35, 40
FTM 791-5422	26			IP 123	10	IP 386	41
FTM 791-5431	34			IP 125	26, 28	IP 388	17, 27
FTM 791-5501	21			IP 131	15	IP 390-A	22
FTM 791-5702	26			IP 132	32	IP 394	16
FTM 791-6503	31			IP 133	32	IP 395	14
FTM 791-6514	31			IP 135	28	IP 398	18
FTM 791-8001	4			IP 136-1	32	IP 409	16
				IP 137	25	IP 410	14
				IP 138	17	IP 411	14
				IP 139	30	IP 412	23
				IP 142	17, 33	IP 421	29
				IP 143	21	IP 423	19
				IP 146	29	IP 440	21
				IP 154	18	IP 441	25
				IP 156	15	IP 459-1	37
				IP 157	27	IP 460-2-3	37
				IP 158	34	IP 481	16
				IP 160	5, 6	IP 491	11, 12
				IP 161	14	IP 492	11, 12
				IP 163	26	IP 506	37
				IP 170	11, 12	IP 515	36
				IP 177	30	IP 516	36
				IP 179	35, 40	IP 520	36
				IP 182	21	IP 523	12
				IP 183	29	IP 524	12
				IP 189	5	IP 540	15
				IP 190	5, 6		
				IP 195	10		



German Standards

<u>DIN</u>	<u>page</u>
DIN 1995 (obs.)	5, 6, 36
DIN 1996-5	38, 42
DIN 1996-15	36
DIN 12 685-2	21
DIN 12 804	14
DIN 51 346	25
DIN 51 347	31
DIN 51 348	29
DIN 51 350-1-6	31
DIN 51 351	24
DIN 51 352-1+2	27, 28
DIN 51 352-2	26
DIN 51 354	31
DIN 51 360-1+2	26, 28
DIN 51 365	24
DIN 51 366	8
DIN 51 367	25
DIN 51 368	25
DIN 51 369	26
DIN 51 372 (obs.)	8
DIN 51 373	27, 28
DIN 51 382	25
DIN 51 385	25
DIN 51 394 (obs.)	27
DIN 51 398	25
DIN 51 400-2-3	20
DIN 51 400-8 (obs.)	20
DIN 51 403-2	17
DIN 51 406 (obs.)	18
DIN 51 408-1	20
DIN 51 411	23, 24
DIN 51 421 (obs.)	16
DIN 51 423-2	42
DIN 51 428-1 (obs.)	16
DIN 51 432	42
DIN 51 538	27
DIN 51 550 (obs.)	8
DIN 51 551-1	18
DIN 51 554-1+2	26
DIN 51 558-1+2	30
DIN 51 559-2	32
DIN 51 560 (obs.)	38
DIN 51 561 (obs.)	7
DIN 51 562-1+2	7
DIN 51 568	24
DIN 51 569 (obs.)	7
DIN 51 571-2 (obs.)	34
DIN 51 575	26
DIN 51 576	21
DIN 51 578 (obs.)	23, 24
DIN 51 579	35, 40
DIN 51 580	35
DIN 51 583 (obs.)	24
DIN 51 587	27
DIN 51 593 (obs.)	24
DIN 51 595	21
DIN 51 597 (obs.)	25
DIN 51 599 (obs.)	23
DIN 51 724-1	20
DIN 51 727	20
DIN 51 750-1-3	4
DIN 51 755-1-2	11
DIN 51 757	5
DIN 51 757-2	6
DIN 51 757-5	14
DIN 51 759-1 (obs.)	18

<u>DIN</u>	<u>page</u>
DIN 51 769-1 (obs.)	21
DIN 51 775 (obs.)	15
DIN 51 780 (obs.)	17
DIN 51 784	15
DIN 51 787 (obs.)	15
DIN 51 793 (obs.)	22
DIN 51 794	12
DIN 51 799 (obs.)	17
DIN 51 801-1	32
DIN 51 802	34
DIN 51 804-2 (obs.)	35
DIN 51 805	33
DIN 51 807-1	34
DIN 51 808	17, 33
DIN 51 809-1 (obs.)	32
DIN 51 809-2 (obs.)	30
DIN 51 811	18
DIN 51 814	34
DIN 51 817	33
DIN 51 825	32
DIN 51 900-1-3	19
DIN 52 004	5, 6
DIN 52 007-1	7, 8
DIN 52 011 (obs.)	36
DIN 52 012 (obs.)	36
DIN 52 013	36
DIN 52 015 (obs.)	37
DIN 52 024	37
DIN 52 025	36
DIN 52 027	37
DIN 53 015	9
DIN 53 211	9
DIN 53 224	9
DIN 53 242-1-4 (obs.)	4
DIN 53 914 (obs.)	24

French Standards

<u>NF</u>	<u>page</u>
FD M07-008	20
NF M07-003	23, 24
NF M07-021	15
NF M07-022 (obs.)	21
NF M07-024	15
NF M07-028	18
NF M07-036	11
NF M41-008 (obs.)	14
NF M41-011 (obs.)	14
NF M41-012	14
NF T20-036	12
NF T20-037	12
NF T30-070 (obs.)	9
NF T60-100 (obs.)	7, 8
NF T60-104 (obs.)	23, 24
NF T60-105	25
NF T60-116 (obs.)	18
NF T60-117	18
NF T60-132	35
NF T65-001 (obs.)	38
NF T66-003 (obs.)	37
NF T66-008 (obs.)	36
NF T66-009	11, 12
NF T66-011 (obs.)	37
NF T66-026 (obs.)	36
NF T66-047 (obs.)	37
NF T78-102	16

European Standards

<u>EN</u>	<u>page</u>
EN 10 (obs.)	15
EN 58	4
EN 116	16
EN 535 (obs.)	9
EN 924	11, 12
EN 1238	36
EN 1426	38
EN 1427	36
EN 1428	41
EN 1871-F	36
EN 12 591	38
EN 12 593	36
EN 12 595	7, 8
EN 12 596	8
EN 12 606-1	37
EN 12 607-1-3	37
EN 12 662	21
EN 12 846	38
EN 13 016-1+2	16
EN 13 179	36
EN 13 303	37
EN 13 358	37
EN 13 398	36
EN 13 589	36
EN 13 703	36
EN 13 880-2	40
EN 14 112	17, 41
EN 14 214	16
EN 14 370	24
EN 14 522	12
EN 14 771	38
EN 15 751	17
EN 16 091	17
EN 22 719 (obs.)	11, 12, 13
EN 23 015	25
EN 24 260	20
EN 26 246 (obs.)	15
EN 60 156	29
EN 61 000-3-2	41
EN 61 000-3-3	41
EN 61 125-A-B-C	27



International Standards

ISO	page	ISO	page
ISO 1409	24	ISO 10 307-1+2	22
ISO 1516	11, 12	ISO 10 337	41
ISO 1523	11, 12	ISO 10 370	18
ISO 1928	19	ISO 11 007	34
ISO 2049	23, 24	ISO 12 205	17, 27
ISO 2083 (obs.)	21	ISO 12 937	41
ISO 2137	35, 40	ISO 13 736	11, 12
ISO 2160	14, 18	ISO 13 758	14
ISO 2176	32	ISO 15 267	11, 12, 13
ISO 2207	32	ISO 15 528	4
ISO 2431	9	ISO 20 623	31
ISO 2592	11, 13		
ISO 2719-A+B	11, 12, 13		
ISO 2908 (obs.)	34		
ISO 2977-A+B+D+E	15		
ISO 3007	16		
ISO 3007-A-B-D	15		
ISO 3012	20		
ISO 3013	16		
ISO 3014	18		
ISO 3015	25		
ISO 3016	25		
ISO 3104	7, 8		
ISO 3105	7, 8		
ISO 3170	4		
ISO 3405	10		
ISO 3675	5		
ISO 3679	11, 12		
ISO 3680	11, 12		
ISO 3733	41		
ISO 3734	22		
ISO 3735	22		
ISO 3771	30		
ISO 3837	15		
ISO 3838	5, 6		
ISO 3841	32		
ISO 3987	26		
ISO 3993	14		
ISO 4256	14		
ISO 4260	20		
ISO 4262	18		
ISO 4263	17		
ISO 4625	36		
ISO 5275	21		
ISO 5661	42		
ISO 6245	26		
ISO 6246	15		
ISO 6247	29		
ISO 6249	19		
ISO 6251	14		
ISO 6293	32		
ISO 6295	24		
ISO 6299	32		
ISO 6614	23		
ISO 6615	18		
ISO 6617-1+2	27, 28		
ISO 6886	41		
ISO 6889	24		
ISO 7120	28		
ISO 7536	17		
ISO 7624	27		
ISO 8819	14		
ISO 9001	8		
ISO 9029	41		
ISO 9038	12		
ISO 9120	29		

Japanese Standards


JIS	page
JIS C2101	4, 26, 27, 29, 30, 42
JIS K2207	5, 7-9, 23, 36-38, 40
JIS K2208	38
JIS K2220	18, 29, 31, 33, 35
JIS K2234	29
JIS K2235	32, 34, 35
JIS K2240	14, 20
JIS K2241	24
JIS K2246	28
JIS K2249	5, 6
JIS K2254	4, 10
JIS K2256	15
JIS K2258	14-16
JIS K2261	15
JIS K2265	5, 11-13
JIS K2269	25
JIS K2270	18
JIS K2272	26
JIS K2275	41
JIS K2276	16, 17, 20, 21
JIS K2279	19
JIS K2283	7, 8
JIS K2287	17
JIS K2288	16
JIS K2501	30
JIS K2503	22, 32
JIS K2510	28
JIS K2513	18
JIS K2514	27, 28
JIS K2518	29
JIS K2519	31
JIS K2520	23
JIS K2536	15
JIS K2537	18
JIS K2541	20
JIS K2541-1	20
JIS K2580	23, 24
JIS K2601	21, 22, 25

International Standards

	page
API 2542	22
API 2548	22
CEC L-07-A-71	31
CEC L-14-A-95	25
CEC L-36-A-90	25
CEC L-40-A-93	29
CEC L-45-T-99	31
CEI 79-4	12
CFR 49-173.120	12
DAB 10-V.5.8.1	35, 40
DIS 3658 (obs.)	6
IEC 156 (IEC 60156)	29
MIL-G-10924	35
NAC TM0172-2001	28
PV 1417	31
PV 1437	31
PV 1444	31
SHELL SMS 658	35, 40
SHRP 1001	38
SHRP 1007	38
SHRP B-003	38
SIS 155 122 (obs.)	16
Superpave®	11, 13, 37, 38
UN Class 3 L.3	12
UN Class 3 N	12
USEPA SW-846-9079	30
VDE 0370-5	29
VDI 2263	13
VV-L-791e-531	33
79/831/EWG	12, 13

Abbreviations

AASHTO	American Association of State Highway & Transportation Officials
API	American Petroleum Institute
ASTM	American Society for Testing & Materials
BS	British Standards Institution
CEC	Coordinating European Council (for Transportation Fluids Performance)
DIN	Deutsches Institut für Normung
EN	Comité Européen de Normalisation
FTM	Federal Test Method (US)
IEC	International Electrotechnical Commission
IP	Institute of Petroleum Standards
ISO	International Standardization Organisation
JIS	Japanese Industrial Standards Committee
MIL	Military Standards (US)
NF	Institut Française de Normalisation
SAE	Society of Automotive Engineers (US)
SHRP	Strategic Highway Research Program (US)



Certificate

Standard **ISO 9001:2008**

Certificate Registr. No. 01 100 005115

TÜV Rheinland Cert GmbH certifies:

Certificate Holder:



Petrotest Instruments GmbH & Co KG
Ludwig-Erhard-Ring 13
D - 15827 Dahlewitz

Scope:

design, manufacture, distribution and repair of petroleum testing equipment and texture analysers

An audit was performed, Report No. 005115. Proof has been furnished that the requirements according to ISO 9001:2008 are fulfilled.

The due date for all future audits is 19-09 (dd.mm).

Validity:

The certificate is valid from 2009-10-09 until 2012-09-19.
First certification 1994

Cologne, 2009-10-09

S. Gebel
TÜV Rheinland Cert GmbH *)
Am Grauen Stein · 51105 Köln



TGA-ZM-58-95-00

10093 4.08 E A4 © TÜV, TÜEV and TUV are registered trademarks. Utilisation and application requires prior approval.

www.tuv.com





You need more details?

The complete product range will be found on our **HOME PAGE** www.petrotest.com where it will be updated immediately.

Search - Functions:

Search by Title, Standard, Order Number, or Product- & User Group.




Leaflets:

Available leaflets as PDF-download files.



Occasions:

Are you looking for a bargain? You may find a previous model or a demo-unit here.




A complete set of the detailed datasheets is available on CD, as a PDF-catalog.

Our catalog includes a comfortable search-function and PDF-files of our leaflets.

MANUFACTURER OF PETROLEUM TEST EQUIPMENT SINCE 1873



PETROTEST BUILDING



DEVELOPMENT

ISO 9001
CERTIFIED
PROCESSES...

...CONSOLIDATED
AT PETROTEST®!




PRODUCTION



QUALITY CONTROL



TRAINING



STORAGE

Presented by:

We reserve the right of changes without prior notice

petrotest®

98-1100 10/09

Ludwig-Erhard-Ring 13 · 15827 Dahlewitz - Germany · Tel.: +49 (0) 33708 / 56-300 · Fax: +49 (0) 33708 / 56-556
E-Mail: info@petrotest.com · Website: www.petrotest.com